# Precise

## Dosing and Dispensing

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

### Operating Manual

Page 2

Peristaltic Pumps Hei-FLOW Value, Hei-FLOW Advantage, Hei-FLOW Precision

Operating manual must be read before initial start-up.

Please follow the safety instructions provided.

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<b>About this Document</b>	<b>4</b>
Version and variants .....	4
About this manual .....	4
Symbols and keywords.....	5
<b>Safety Instructions</b>	<b>6</b>
General safety instructions.....	6
EU Declaration of Conformity.....	6
Intended use.....	6
Installation / electrical safety .....	7
Qualifications of personnel .....	7
Obligations of the operator .....	7
Installation site .....	7
Modifications to the device.....	8
Safety of personnel.....	8
Safety during use .....	9
Disposal .....	9
<b>Device Description</b>	<b>10</b>
Device overview .....	10
Hei-FLOW Value: overall view .....	10
Hei-FLOW Advantage: overall view.....	10
Hei-FLOW Precision: overall view .....	11
Hei-FLOW Precision: control panel.....	11
<b>Start-Up</b>	<b>12</b>
Setting up device .....	12
Switch main switch on and off.....	12
Switch unit on .....	12
Switch unit off.....	12
Connect/disconnect power cord .....	13
<b>Operation</b>	<b>14</b>
Starting, stopping pump.....	14
Hei-FLOW Value.....	14
Hei-FLOW Advantage, Hei-FLOW Precision.....	14
Change direction of rotation.....	16
Additional functions Hei-FLOW Precision .....	17
Set pump head.....	18
Set tube .....	18
Set flow.....	19
Set volume .....	19
Calibrate pump.....	20
Set intermittent dosing.....	21
Set back stroke .....	22

Interfaces.....	22
Analog Interface .....	23
Digital Interface .....	24
Digital interface commands.....	25
<b>Cleaning and Maintenance</b>	<b>27</b>
Cleaning.....	27
Maintenance .....	27
<b>Troubleshooting</b>	<b>28</b>
Possible errors and how to resolve .....	28
<b>Assembly</b>	<b>29</b>
Electrical connection .....	29
Mount pump heads .....	30
Single-channel pump heads SP standard and SP vario .....	30
Single-channel pump head SP quick .....	31
Multi-channel pump heads .....	32
Mount tubes: Single-channel pump heads.....	34
Tubes: pump head SP standard .....	34
Tubes: pump head SP vario.....	35
Tubes: pump head SP quick .....	37
Mount tubes: Multi-channel pump heads.....	38
Cassettes and tubes: „Cassette small“ .....	39
Cassettes and tubes: „Cassette medium“ .....	41
Cassettes and tubes: „Cassette large“ .....	42
Mount several pump .....	43
<b>Disassembly and Storage</b>	<b>43</b>
Dismantling, transportation and storage .....	43
Dismantling .....	43
Transportation and storage .....	43
<b>Accessories and Spare Parts</b>	<b>44</b>
Scope of delivery .....	44
Accessories.....	44
<b>Attachments</b>	<b>45</b>
Technical data .....	45
Table tubes .....	46
Performance range pump heads.....	48
<b>Service</b>	<b>49</b>
Contact / Technical Service .....	49
Warranty .....	50
Confirmation of condition.....	51

### Copyright

This operating manual is protected by copyright and is exclusively intended for the internal use of the purchaser.

Handing over this operating manual to third parties, reproducing it in any type or form –in part or in full – and utilizing and/or communicating its contents is prohibited without prior written approval from the manufacturer.

Any violations will result in an obligation to provide compensation for damages. Other claims remain unaffected.

## ➤ Version and variants

### Version

This manual describes the function, operation and maintenance of the Hei-FLOW Value, Hei-FLOW Advantage and Hei-FLOW Precision peristaltic pumps.

Version	Alteration date
1.2	10/2018

### Variants

The device is available in different variants. Certain features or functions are only available in specific product variants. These variants are described in this operating manual.

## ➤ About this manual

This operating manual is an integral part of the device described here.

- Carefully read through this operating manual and observe all of the safety instructions and warnings notices.
- Ensure that every user has carefully read the operating manual before using the device for the first time.
- Store this operating manual so that it is accessible to all users at all times.
- Pass the operating manual on to any subsequent owner.

The current version of this operating manual can also be found at:

[https://heidolph-instruments.com/en/service/downloads/operation-manuals#Peristaltic Pumps \(Hei-FLOW Value-Advantage-Precision\)](https://heidolph-instruments.com/en/service/downloads/operation-manuals#Peristaltic%20Pumps%20(Hei-FLOW%20Value-Advantage-Precision))

In addition, also observe the local regulations at the installation site.

## ➤ Symbols and keywords

Standardized symbols and signal words are used in this operating manual to warn against any dangers and provide important instructions. These instructions must be strictly observed to avoid accidents and damage.

The following terms and basic symbols are used:

Symbol	Additional signal words / explanation
<b>Warning signs</b> 	The yellow triangle indicates hazardous situations. It is used in combination with the following signal words: <b>DANGER:</b> Indicates a hazardous situation which, if not avoided, will result in serious injury or death. <b>WARNING:</b> Indicates a hazardous situation which, if not avoided, may result in serious injury or death. <b>CAUTION:</b> Indicates a hazardous situation which, if not avoided, may result in property damage and minor or moderate injury.
<b>Prohibitory signs</b> 	Prohibited: The red circle indicates a situation that should be avoided under all circumstances and which, if not avoided, may result in serious injury or death.
<b>Mandatory signs</b> 	Must be observed: The blue circle indicates important information. Please observe this information to avoid any property damage.

### Other symbols used:

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a.	List of variants
b.	
c.	

## ➤ General safety instructions

### EU Declaration of Conformity

This device complies with the following EC-Directives\*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

\*Also see the attachment "EU Declaration of Conformity".

The device has been manufactured according to state-of-the-art technology and in compliance with recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure that the operating manual is available at all times.

#### The device may only be used under the following conditions:

- Only operate the device if it is in a technically perfect condition.
- Only operate the device if you are properly aware of the risks and required safety measures.
- Only operate the device in accordance with the instructions given in this operating manual.
- If any information is ambiguous or missing, ask your superior or contact the manufacturer.
- Do not operate anything on the device without authorization.
- Only use the device in accordance with its intended use.

### Intended use

The device may only be operated by authorized personnel.

The device is suitable for the following use:	Hei-FLOW Value Hei-FLOW Advantage	Dispensing of fluids
	Hei-FLOW Precision	Dispensing and dosing of fluids
Operation of the device is suitable in one of the following laboratories:	<ul style="list-style-type: none"> <li>▪ Chemistry</li> <li>▪ Pharmacy</li> <li>▪ Biology</li> <li>▪ Environmental analytics</li> <li>▪ Basic research</li> <li>▪ Any other research laboratory</li> </ul>	

#### Reasonable foreseeable misuse

Any use which deviates from the device's intended use is considered to be improper. The manufacturer is not liable for damage that occurs as a result. The risk is borne by the operator alone.

### Installation / electrical safety

- The device may only be connected up if the available power supply voltage corresponds to the information stated on the rating plate for the device.
- The power supply connection must be easily accessible at all times.
- An RCD circuit breaker (residual current circuit breaker) must be fitted to the electrical system in the building so that the device is disconnected from the power circuit in the event of a fault.
- Repairs may only be performed by electricians authorised by Heidolph Instruments.
- Never operate the device with a damaged power cord.
- Always turn the device OFF and disconnect the power cord before carrying out any maintenance or repair work.

### Qualifications of personnel

- The device may only be operated by trained personnel.
- The device may only be operated by persons who have been instructed and supervised in its proper use by trained specialist personnel.
- The device may only be operated by specialist personnel who are above the legal minimum age.
- Other persons may only work on the device under the constant supervision of experienced and trained specialist personnel.
- This operating manual must be read and understood by all persons working with the device.
- Personnel must receive safety training that ensures responsible and safe working practices.

### ➤ Obligations of the operator

#### Installation site

- The device must be installed in a suitable location.
- The device must be installed on a firm and stable surface.
- Ensure that the device and all of its components are easily accessible at all times.
- Maintain a sufficient safety clearance to, in particular, any moving and/or hot equipment components.
- It is not permitted to place or store any objects such as accessories, tools or chemicals within this clearance area during operation.
- All screw connections must be securely tightened.
- It is not permitted to operate the device near to highly flammable or explosive substances.
- Operate the device in conjunction with an extractor hood, if working with potential harmful media (see DIN EN 14175 and DIN 12924).

## Safety Instructions

- The device is designed for indoor use only and under the following ambient conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity
	32 - 40 °C decreasing linearly to a maximum 50 % relative humidity
Installation altitude	0 - 2,000 m above sea level
Contamination level	2
Overvoltage category	II
Permissible supply deviations	± 10 %

- If the device is operated in corrosive atmospheres, the service life of the device will decrease based on the concentration, duration and frequency of the exposure to the corrosive atmosphere e.g. concentrated hydrochloric acid (HCl).

## Modifications to the device

- It is not permitted to make any unauthorised modifications or changes to the device.
- Do not attach or install any parts that have not been approved by the manufacturer.
- Unauthorized modifications or changes will void the EC Declaration of Conformity for the device and operation of the device will no longer be permitted.
- The manufacturer is not liable for any damage, dangers or injuries that result from unauthorized modifications and changes or due to the non-observance of the instructions in this manual.

## Safety of personnel

- Ensure that the device is only operated by qualified specialist personnel and trained employees.
- Observe the following instructions to avoid any personal injuries and property damage:
  - Laboratory regulations
  - Accident prevention regulations
  - Hazardous Substance Act
  - Other generally accepted occupational health and safety regulations
  - Local regulations

### ➤ Safety during use

- Wear the appropriate clothing when working on the device (protective clothing, protective glasses and, if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not carry out work with naked flames in the vicinity of the device (risk of explosion).
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range ( $3 \cdot 10^{11}$  Hz to  $3 \cdot 10^{15}$  Hz).
- Do not operate or assemble appliances in the vicinity of the device which constitute emission or radiation sources for ionizing radiation or in the ultrasonic range.
- Do not operate the device where adiabatic compression or shock waves might occur (shock wave combustion).
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediatly.
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediatly.
- Eliminate errors immediately.
- Always switch the device OFF after use.
- Avoid putting pressure on the display when you are not operating the device.
- Do not use abrasive material to clean the surface. Only wipe with damp cloth.

### ➤ Disposal


- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device in accordance with the relevant national regulations.
- Dispose of the packaging material in accordance with the appropriate national regulations.


## ➤ Device overview

### Hei-FLOW Value: overall view

### Hei-FLOW Advantage: overall view


# Device Description

## Hei-FLOW Precision: overall view

## Hei-FLOW Precision: control panel

### ➤ Setting up device

**Warning:**  
**Risk of slipping device!**

In operation the device might slide off the table top.

The surface must be smooth and clean and should be properly leveled.

- ➔ Locate the unit on a stable, horizontal surface.
- ➔ Clean surface and feet with a damp cloth and ethanol regularly.

### ➤ Switch main switch on and off

The ON/OFF switch is situated on the back side of the unit.

- ✓ ON/OFF switch is in „power off“ position.

### Switch unit on

- ➔ Press ON/OFF switch.
- ✓ ON/OFF switch is in „power on“ position.

**Hei-FLOW Advantage**

- ✓ The indicator light of actual direction illuminates.

**Hei-FLOW Precision**

- ✓ The settings are displayed in following order:
  - set pump head (e. g. „Pu 1“)
  - gear ratio (e. g. „5“)
  - set inner tubing diameter (e. g. „3,1“)
  - possible back hub (z. B. „r on“)
- ✓ The indicator light of actual direction illuminates.
- ✓ Unit is switched on.

### Switch unit off

- ➔ Press ON/OFF switch.
- ✓ ON/OFF switch is in „power off“ position.
- Hei-FLOW Advantage, Hei-FLOW Precision**
  - ✓ All indicator lights power off.
- ✓ Unit is switched off.

### ➤ Connect/disconnect power cord

The three-wire socket is firmly connected to the unit.

The connection is at the rear side of the unit.

- ✓ The ON/OFF switch is in the "Off" position.
- Connect the power cord with a grounded socket.
- Disconnect the plug if the unit is not in use.

## Starting, stopping pump

### Hei-FLOW Value

#### Set speed

Speed is set via rotation knob.

- ✓ Unit is switched on.
- ✓ Required direction of rotation is activated.
- Turn knob to set speed.


#### Start rotation

- Press rotation switch on the right side.
- ✓ Indicator light of rotation switch illuminates.
- ✓ Pumping starts.

#### Stop rotation

- Press rotation switch on the left side.
- ✓ Indicator light of rotation switch powers off.
- ✓ Pumping stops

Rotation switch

Rotation on

Rotation off

A power cut will stop the pump. When power is returned the indicator light of the rotation switch will start blinking.

Rotation must be restarted by switching the rotation switch off and on again.

## Hei-FLOW Advantage, Hei-FLOW Precision

#### Caution: Possible damage of device

If device sockets are not covered, fluids can get in contact with Interface connectors. Penetrating fluids can lead to malfunctions of interface connection and could initiate a short circuit.

Sockets at the backside of device must be covered at all times. Close sockets with attached sealing cap or connect interface cable.

## Operation

### Set speed: Hei-FLOW Advantage

Speed is set via rotation knob.

- ✓ Unit is switched on.
- ✓ Required direction of rotation is activated.
- Turn knob to set required speed.
- Switch rotation on (see underneath).


### Set speed: Hei-FLOW Precision

- Hei-FLOW Precision 01: 5 - 120 rpm
- Hei-FLOW Precision 06: 24 - 600 rpm

Speed is set via selection key „Speed rpm”.

- ✓ Unit is switched on.
- ✓ Required direction of rotation is activated.
- Press selection key repeatedly, till indicator light „Speed rpm” illuminates.
- Increase speed via arrow key up .
- Decrease speed via arrow key down .
- ✓ Set speed will be displayed as rpm.


### Switch rotation on

- Press rotation key.
- ✓ The indicator light of rotation key will illuminate.
- ✓ Pumping starts.

### Switch rotation off

- Press rotation key.
- ✓ The indicator light of rotation key will power off.
- ✓ Pumping stops.

### Rotation key

Feed rate depends on media, inner diameter and material of tubing and speed.

To determine speed see chapter „Technical data, Performance range”.

## Operation

### Charge and discharge tubes

With the key „max“ charging and discharging of the tubes can be sped up without changing any settings.

- ✓ Unit is switched on.
- ✓ Rotation is switched off.
- ✓ Required direction of rotation is activated.
- Press key „max“ and hold pressed.
- ✓ Pump is running with maximum speed.

#### Hei-FLOW Precision:

- ✓ The display shows „Full“.
- ✓ Release key „max“.
- ✓ Pumping stops.

A power cut will stop the pump. All settings will be saved.  
With power return the device is switched on automatically. Rotation must be started manually.


## ➤ Change direction of rotation

Direction of rotation can only be changed if rotation is switched off.  
A warning signal will be audible when trying to change direction during operation.

### Change direction: Hei-FLOW Value

- ✓ Rotation is switched off.
- Press direction switch on the right side.
- ✓ Pump will operate clockwise.
- Press direction switch on the left side.
- ✓ Pump will operate counter-clockwise.


## Operation

### Change direction: Hei-FLOW Advantage, Hei-FLOW Precision

- ✓ Rotation is switched off.
- ✓ The left control light of direction key is illuminated.
- ✓ Rotation is switched off.
- Press direction key.
- ✓ The right indicator light of direction key illuminates.
- ✓ The pump will operate clockwise.
- Press direction key.
- ✓ The left indicator light of direction key illuminates.
- ✓ The pump will operate counter-clockwise.

Counter-clockwise

Clockwise

### ➤ Additional functions Hei-FLOW Precision

On top of the standard functions you can set parameters that allow consistent process quality.

Following parameters can be set:

- „Speed rpm“  
sets the speed of pumping;  
flow per minute “Flow ml/min” will be adjusted automatically
- „Tube mm” (inner diameter)  
allows minimum and maximum flow per minute  
must be set correctly for any other parameter to be adjusted correctly
- „Flow ml/min”  
sets the flow per minute;  
speed “Speed rpm” will be adjusted automatically
- „Vol ml”  
sets the volume to be pumped in total;  
pump will be stopped automatically when volume is reached
- „Pause sec/min/h”  
sets the time per pause:  
intermittent dosing allows a permanent repetition of dispensing a defined volume

Correct operation of all functions requires the correct setting of pump head and inner tube diameter in use.

When purchasing a package with pump head, the pump head is set correctly by the factory.


## Set pump head

Factory pre-setting is SP quick (Pu 1).

Display	Pump head		Hei-FLOW Precision 01	Hei-FLOW Precision 06
Pu 0	Multi-channel	C 8	✓	-
Pu 1	Single-channel	SP quick	✓	✓
Pu 2	Single-channel	SP standard / SP vario	✓	✓
Pu 3	not assigned	-	✓	-
Pu 4	not assigned	-	✓	-
Pu 8	Multi-channel	C 12	✓	-
Pu 9	Multi-channel	C 4	✓	-

- Switch unit off at main switch.
- Press „Cal“ key and hold pressed.
- Switch unit on at main switch.
- ✓ Display shows actual setting of pump head.
- Release „Cal“ key.
- Select required pump head via arrow keys  or .
- Confirm selection by pressing „Select“ key.

## Set tube

Under "Tube mm" set the inner diameter of the tube. Possible settings vary depending on the pump head and the corresponding tube in use. (For tubes and settings please see chapter „Technical data, Selection table tubes“.)

- ✓ Unit is switched on.
- Press „Select“ key repeatedly, till indicator light „Tube mm“ illuminates.
- Select required inner diameter of tube in use via arrow keys  or .
- Confirm selection by pressing „Select“ key.

A blinking indicator light „Vol ml“ during setting of tube diameter shows that the device adjusts volume to required minimum of new diameter.

## Set flow

By setting the flow you define how many ml per minute should be dispensed.

Depending on the tube size different ranges of flow.  
(Please refer to chapter "Technical data, Table tubes".)

When using „Cassette small“ with small tubes and inner Ø from 0,2 mm - 0,9 mm volume will be displayed as **µl** instead of **ml**.

- ✓ Unit is switched on.
- ✓ Pump head and inner diameter of tube is set correctly.

- Press „Select“ key repeatedly, till indicator light „Flow ml/min“ illuminates.
- ✓ Display shows actual setting of flow ml per minute.
- Select required flow ml per minute via arrow keys  or .
- Confirm selection by pressing „Select“ key.
- ✓ The corresponding speed will be adjusted automatically.

## Set volume

By setting the volume you define how many milliliters in total should be dispensed.

Depending on the tube size used the volume can be set from 0,1 to 9,999 ml.

When using „Cassette small“ and tubes with an inner Ø between 0,2 mm and 0,9 mm the volume will be displayed as **µl** instead of **ml**.

Model Hei-FLOW Precision 06 is limited to a maximum speed of 300 rpm.

- ✓ Unit is switched on.
- ✓ Pump head and inner diameter of tube is set correctly.

- Press „Select“ key repeatedly till indicator light „Vol ml“ illuminates.
- Select required volume ml via arrow keys  or .
- Confirm selection by pressing „Select“ key.
- ✓ Pump will automatically stop when volume is reached.

Speed and tube diameter result in a minimal volume. Lower volume can only be performed by using smaller tubes.

Large tubes and too little flow will result in an error feedback (see chapter "Troubleshooting, Possible errors and how to resolve").

## Calibrate pump

Tubes have individual manufacturing tolerances. If high dosing accuracy is required the pump has to be calibrated for the tube in use by flow and/or by volume.

- ✓ Unit is switched on.
- ✓ Media to be dispensed is connected to pump.
- ✓ Pump head and inner diameter of tube are set correctly.

### Calibrate pump with flow

- ✓ Indicator light „Flow ml/min” is illuminated.
- ✓ Required flow ml per minute is set.
- Press „Cal” key.
- ✓ Display starts blinking.
- Press rotation key.
- ✓ The pump dispenses for exactly 1 minute.
- Measure dispensed volume.
- Set measured flow via arrow keys  or .
- Press „Cal” key for a minimum of 3 seconds to confirm.
- ✓ Pump is calibrated for flow ml per minute.

Verify set flow after calibrating. If applicable adjust flow.

### Calibrate pump with volume

- ✓ Indicator light „Vol ml” is illuminated.
- ✓ Required volume is set.
- ✓ Speed is set via setting flow ml/min.
- Press rotation key.
- Measure dispensed volume.
- Press „Cal” key.
- Set measured volume via arrow keys  or .
- Press „Cal” key for a minimum of 3 seconds to confirm.
- ✓ Pump is calibrated for volume.

If the calibrated values show a deviation > 30% most probably pump head or tube diameter are not set correctly. Or the tube mounting needs to be adjusted.

In this case cancel calibration via „Select” key. Verify your settings.

## Set intermittent dosing

Pauses can be set from 0.1 seconds to 750 hours.

Time unit	Display	Explanation
Seconds	0.1 - 599.9	600 sec. = 10 minutes
Minutes	10 - 599*	600 min. = 10 hours
Hours	10 - 750h	750 h. = 31.25 days

- ✓ Unit is switched on.
- ✓ Pump head and inner diameter of tube are set correctly.
- ✓ Volume is set.
- Press the „Select“ key repeatedly till the indicator light „Pause sec/min/h“ illuminates.
- Select the required break time between pumping motion via arrow keys  or .
- Press the „Select“ key repeatedly till indicator lights „Vol ml“ and „Pause sec/min/h“ illuminate simultaneously.
- ✓ Intermittent dosing is set for all following processes.
- ✓ The display shows alternating dispensed volume or running break time. The corresponding indicator light blinks.
- ✓ As soon as break time is over the pump starts moving.
- ✓ As soon as set volume is dispensed the pump will stop.


### Set back stroke

To avoid dripping when dosing a short backstroke can be programmed.

- Switch unit off at main switch.
- Press „Cal“ key and hold pressed.
- Switch unit on at main switch.
- ✓ Display shows actual setting of pump head.
- Release „Cal“ key.
- Press „Select“ key repeatedly till Display shows „r on“ or „r off“ .
- ✓ „r on“ switches back stroke on, „r off“, switches back stroke off.
- Confirm selection by pressing „Select“ key.

### ➤ Interfaces

#### **Warning:** **Danger of electric shock!**

If the interface inputs and outputs are supplied with excessively high voltage and are not sufficiently insulated, in the event of a fault, metallic parts (e.g. the housing), will be electrified.

Voltage inputs and outputs above 25 V AC or 60 V DC are deemed safe if separated according to DIN EN 61140, or by double or reinforced insulation in accordance with DIN EN 60730-1 or DIN 60950-1.

Only use shielded connection cables. Connect shield to connector case.

#### **Caution:** **Possible damage of device**

If device sockets are not covered, fluids can get in contact with Interface connectors. Penetrating fluids can lead to malfunctions of interface connection and could initiate a short circuit.

Sockets at the backside of device must be covered at all times.

Use plug connectors protection class IP 67 only. Ensure professional connection.

## Analog Interface

(Hei-FLOW Advantage and Hei-FLOW Precision only)

To operate pump externally an analog interface is located at the backside of device.


Via this interface the speed can be set.

The socket is protected against splash water with a screw cap.

### Pin assignment of analog interface socket

- PIN 1 = +5 V
- PIN 2 = Start/Stop
- PIN 3 = 0-10 V
- PIN 4 = Direction of rotation
- PIN 5 = 0 V (GND)
- PIN 6 = 4-20 mA

Control voltage	Control current	Hei-FLOW Advantage 01 Hei-FLOW Precision 01	Hei-FLOW Advantage 06 Hei-FLOW Precision 06
0 - 0.9 V	0 - 4.9 mA	0 rpm	0 rpm
1 V	5 mA	5 rpm	24 rpm
10 V	20 mA	120 rpm	600 rpm

- TTL-level 0 V = clockwise rotation, TTL-level 5 V = counter-clockwise rotation
- TTL signal 5 V switches rotation on and off. This function is available all the time.
- For triggering direction of rotation and start/stop function for rotation, a 5 volts pin is available in the interface connector which allows potential-free switching.
- The analog interface is recognized at a level of 0.2 volts or a 4 mAmps current.
- **Hei-FLOW Precision:** This mode is confirmed by "EA" in the display. As long as the pump drive is controlled via analog interface, the display is inactive and "Start/Stop" key will be the only active key on the device.


# Digital Interface


**(Hei-FLOW Precision only)**

You can connect the device to a PC.

In addition to an analog interface a RS 232 digital interface is located at the backside of device.


**Digital interface socket**

**Pin assignment**

- PIN 1 — Pin 1
- PIN 2 — Pin 6
- PIN 3 — Pin 2
- PIN 4 — Pin 7
- PIN 5 — Pin 3
- PIN 6 — Pin 8
- Pin 7 — Pin 4
- Pin 8 — Pin 5

**Socket cable to PC**

As long as the RS 232 interface is active, all pump control elements and display will remain active.

## Digital interface commands

Query to pump <sup>1)</sup>	Feedback from pump	Explanation
LEDx?	<ul style="list-style-type: none"> <li>• LED=0000\r\n</li> <li>• LED=0001\r\n</li> </ul>	Select parameter: x=0-7 <ul style="list-style-type: none"> <li>▪ LED off</li> <li>▪ LED on</li> </ul> <hr/> LED0?=Volume ml <hr/> LED1?=Flow ml/min <hr/> LED2?=Tube mm (inner diameter) <hr/> LED3?=Speed rpm <hr/> LED4?=Rotation <hr/> LED5?=Direction of rotation right <hr/> LED6?=Direction of rotation left <hr/> LED7?=Pause sec/min/h (intermittent pause)
DSP?	DSP=xxx\r\n	Display query: show actual value
TAx!	-	Select function key: x=1-7 <hr/> TA1!=Cal (calibrate) <hr/> TA2!=Start/Stop (rotation key) <hr/> TA3!=Change direction of rotation <hr/> TA4!=Max (maximum speed) <hr/> TA5!=Up (arrow key) <hr/> TA6!=Down (arrow key) <hr/> TA7!=Select (select key)
SPx=xxxx!	-	Set intermittent pause: x=0;2;3 <hr/> SP0=0001 - 5999 seconds/10 e.g. 100 milliseconds <hr/> SP2=0010 - 0599 minutes <hr/> SP3=0010 - 0750 hours
SMM=xxxx! <sup>2)</sup>	-	Flow ml respectively µl (tube dependent): xxxx=000.0 - 9999
SDZ=xxxx!	-	Set rotation rpm: xxxx=0000 - 9999 (min/max gear dependent)
SDM=xxxx! <sup>2)</sup>	-	Volume dosing ml e.g. µl (tube dependent): xxxx=000.1 - 9999 (not allowed for calibrating)
SSD=xxxx! <sup>2)</sup>	-	Set inner tube Ø : x=4digit figure xxxx=0000 (see following section "Interface commands: set tube")
-	OK\r\n	Command successful
-	ERROR\r\n	Command failed
-	PO/r/n	Dosing completed

<sup>1)</sup> **Interface parameter RS 232:** 9600 Baud, No Parity, 8 Bit, 1 Stopbit

<sup>2)</sup> Also see chapter „Technical data, Table tube selection“ for tube dependent parameters

- Do not send commands as package, minimum pause 0.1 seconds
- xxxx substitutes a 1 - 4 digit figure

## Interface commands: set tube

<b>Command</b>	<b>Inner tube Ø</b>
SSD=0000!	0.8 mm
SSD=0001!	1.0 mm
SSD=0002!	1.7 mm
SSD=0003!	2.0 mm
SSD=0004!	3.0 mm
SSD=0005!	3.1 mm
SSD=0006!	4.0 mm
SSD=0007!	4.8 mm
SSD=0008!	5.0 mm
SSD=0009!	6.0 mm
SSD=0010!	6.3 mm
SSD=0011!	7.0 mm
SSD=0012!	7.9 mm
SSD=0013!	8.0 mm
SSD=0014!	9.0 mm

### > Cleaning

All surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

**Warning:**  
**Risk of electric shock!**

If the device is still connected to power during cleaning, you might get in contact with live parts due to accidentally entering liquids.

Prior to cleaning:

- ensure that the device is switched OFF
- ensure that the plug is pulled from the electric socket.

**Caution:**  
**Damage to the device**

Accidentally entering liquids will cause the device to fail .

Only wipe with a damp cloth. Never use a soaked cloth.

**Attention:**  
**Risk of damage to the surfaces**

All of the surfaces may be damaged by improper cleaning.

Never use:

- Chlorine bleach or any cleaning agent containing chlorine
- Solvent-bearing agents (e.g. acetone)
- Ammonia
- Abrasive cleaning agents such as scouring agents, cleaning wool or agents with metallic particles

To ensure a stable position clean surface and feet with a damp cloth and ethanol regularly.

### > Maintenance

The unit is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

- Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

#### **Tubes**

Tubes are wear parts. Service life of tubes depends on material, rotational speed and properties of the substance to be fed, e.g. temperature, viscosity, chemical composition, and pressure.

Service life of tubes used in pumps varies between approx. 50 and 500 hours in use.

- Check tubes regularly.
- Replace tubes in the event of visible wear.

#### **Multi-channel pump head**

- Check plastic sleeve over metal bar regularly.
- Replace plastic sleeve in the event of visible wear or latest after approx. 500 hours in use.

## Troubleshooting

### ➤ Possible errors and how to resolve

#### General

Errors	Possible reason	Troubleshooting
Unit will not switch on	No power supply	Check size and shape of plug and compatibility to your electrical socket
Indicators do not illuminate	Fuse defective	Please contact your local Heidolph Instruments distributor
	<ul style="list-style-type: none"> <li>▪ Regulation defective</li> <li>▪ Motor defective</li> <li>▪ Indicator light defective</li> </ul>	Please contact your local Heidolph Instruments distributor
	Tube moves	Adjust tube mounting (see chapter „Assembly, Mount tubes“)
Pump is running but does not dispense	Wrong tube wall thickness	Use adequate tube (see chapter „Assembly, Mount tubes“)
	Gear not working: coupling of pump head not aligned	Push pump head to drive as far as it will go (see chapter „Assembly, Mount pump heads“)
	Tube bed open (pump head SP quick only)	Close tube bed (see chapter „Assembly, Mount tubes“)
	Tube not set correctly (Hei-FLOW Precision only)	Set correct tube diameter (see chapter “Hei-FLOW Precision, additional functions”, section “Set tube”)
	Motor too hot: overheat protection activated	Wait a few minutes and let motor cool down; pump switches on automatically

#### Hei-FLOW Advantage, Hei-FLOW Precision

Error message	Possible reason	Troubleshooting
Indicator light rotation key is blinking	Power cut	Start rotation

#### Hei-FLOW Precision

Error message	Possible reason	Troubleshooting
----	Invalid values while calibrating volume	<ul style="list-style-type: none"> <li>▪ Set rotation rpm according to table (see chapter “Technical data, Performance range pump heads”)</li> <li>▪ Calibrate volume anew (see chapter “Technical data, Table tubes”)</li> </ul>

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph distributor immediately.

## Electrical connection

Repairs to the device must **only be carried out by a qualified electrician approved by Heidolph Instruments**. Improper repairs could result in serious dangers. Please contact your specialist dealer if the device needs to be repaired.

### Warning: Risk of electric shock

The device must be connected to an earthed power socket.

The rated voltage for the device and the supply voltage must be the same. The rated voltage for the device can be found on the rating plate on the rear of the device.

The device is earthed via the supplied power cord when delivered. If the original plug is replaced, it is essential that the new plug is connected to a protective earth!

Rating plate

ON/OFF  
switch-Schalter

### Colour coding for the electrical connection:

#### Europe

**GREEN/  
YELLOW**

PE: Protective earth  
**(earth)**

**BLUE**

N: Neutral conductor

**BROWN**

P: Phase

#### North America

**GREEN**

PE: Protective earth  
**(earth)**

**WHITE**

N: Neutral conductor

**BLACK**

P: Phase

The device is connected to the power supply with the supplied power cord. The device socket is located on the rear of the device.

If you want to use the device in a country with a different type of plug:

- The supplied plug may only be replaced by a qualified electrician.
- If you want to use an adapter, it must be earthed and approved in accordance with local regulations.
- The cable, adapter and plug must correspond at least to the power stated on the rating plate.

## Mount pump heads

**Warning:**  
**Risk of finger crushing!**


When mounting pump heads you can crush your fingers if rotation is switched on accidentally.

Always switch device off via main switch while mounting pump heads.


### Single-channel pump heads SP standard and SP vario

Components:

- 1x Pump head
- 1x Flange
- 1x Guide aid
- 2x Knurled screws
- 2x Short tube sections (adaptors for tubes  $\varnothing$  0,8 mm and  $\varnothing$  1,7 mm)
- 1x Allen key

- Insert guide aid into flange.
- Position flange with guide aid in front of pump head in a way that guide will be pushed into coupling.
- ✓ Allen screw shows upwards.
- Tighten flange on both sides with knurled screws onto pump head.
- Remove guide aid.


- Remove protection shield from front of pump head.

**Warning:**  
**Risk of finger crushing!**

When pump head is not covered, you can crush your fingers on rotating parts.


Always mount protection shield before switching on and operating device.

## Assembly

- Push pump head with coupling heading into flange, till flange and pump head are aligned.
- Eventually turn rotor till couplings interlock.
- Position pump head so it fits to your set-up of vessels and fixate with allen screws.
- Mount protective shield.

### Single-channel pump head SP quick

- Move clamping lever from right to left.
- ✓ Tube bed of pump head is open.

Clamping lever

Tube bed open

- Insert coupling of pump head into coupling of pump drive.
- Eventually turn rotor till couplings interlock.

Coupling pump drive

Coupling pump head

- Hand-tighten pump head with the 2 knurled screws (tighten up with a screwdriver or coin if necessary).

Knurled screws


# Multi-channel pump heads

## Components

### Adaptor:

- 1x Adaptor
- 2x Cross headed screws

Adaptor

Front side      Back side      Screws

### Pump head C 8

- 1x Pump head
- 1x Spare plastic sleeve for metal bar
- 4x Allen screws
- 1x Allen key
- 2x Spare feet

Pump head C 8      Plastic sleeve      Allen key

Spare feet      Allen screw


### Pump head C 4

- 1x Pump head
- 4x Allen screws
- 1x Allen key

Pump head C 12      Allen key

### Pump head C 12

- 1x Pump head
- 4x Allen screws
- 1x Allen key
- 2x Spare feet

Spare feet      Allen screw

## Assembly

- Insert coupling of adaptor into coupling of pump drive.
- Insert supplied screws through holes of adaptor.
- Tighten screws with a cross headed screw driver.

**Caution:**  
**Possible damage of device (Pump head C 8 only)**

Without plastic sleeve mounted on metal bar the pump head can be damaged.


Ensure plastic sleeve is mounted on metal bar for operation.

Pumphead C8 only:  
plastic sleeve


- Align metal bar (with plastic sleeve) to mounting slot of adaptor.
- Push pump head onto adaptor in a way that metal bar glides into mounting slot of adaptor.
- Level pump head via adjusting foot for support.
- Insert allen screws into 4 holes and tighten with supplied allen key onto adaptor.

Holes for allen screws

- Tighten screw nut at adjusting foot of pump head.

### ➤ Mount tubes: Single-channel pump heads

**Warning:**  
**Risk of crushing!**

While mounting and demounting tubes you can crush your fingers, if rotation starts accidentally.

Always switch device OFF at main switch before mounting tubes.

Only use tubes with appropriate wall thickness. Using an inadequate tube can damage tube or pump might not operate.


To extend life-cycle of tubes pre-tension tube after mounting. Pull slightly and simultaneously at both ends of tube.

Pre-tension guarantee the rollers running centered over tube.

You will find corresponding tube sizes to variants of pump heads in chapter „Technical data, Table tubes“ as well as in our general catalogue.

### Tubes: pump head SP standard

- ✓ Protection shield is removed.
- Loosen clamping screw to take off clamp.
- Mount tube in a curve of 180° into tube bed.
- Turn rotor by hand to engage tube between rollers and tread.
- Turn rotor repeatedly till tube is crushed by total width of rollers.
- Secure clamp with clamping screws.
- Mount protection shield and secure with knurled screws.
- Check correct mounting by running a short test.

## Tubes: pump head SP vario

Additional components:

1x Aid pin

Pump head SP vario allows operating tubes of different wall thickness.


- ✓ Remove protection shield.
- Remove clamping screw to remove clamp.

Aid pin

- Turn vario wheel via aid pin to maximum wall thickness.
- Mount tube in a curve of 180° into tube bed.
- Turn rotor by hand till tube lies between roller and tread.
- Turn vario wheel via aid pin to minimum wall thickness till you sense a noticeable resistance.
- Remove aid pin.
- Tighten fixing screw to hold vario wheel in this position.

- Turn rotor repeatedly to ensure tube is crushed by rollers in whole width.
- Attach clamp with clamping screws.
- Mount protection shield and secure with knurled screws.
- Check correct mounting by running a short test.


The right pressure is reached if fluid is dispensed during operation but stops passing through when pump is stopped.

### Small tubes, inner diameter 0.8 mm and 1.7 mm

Tubes with an inner diameter of 0.8 mm and 1.7 mm are delivered with 2 additional tube cuttings.

Push both tube cuttings onto the dispenser tube.

- Place tube in pump head with cuttings positioned in clamping sleeve.
- ✓ The cuttings ensure save clamping of the small tube.



### Tubes: pump head SP quick

Pump head SP quick features a quick-clamp mechanic for clamping tubes.

- Move clamping lever to the left to open tube bed.
- Turn both clamping screws counter-clockwise as far as they will go for maximum tube size.

- Place tube into open tube bed.
- Move clamping lever cautiously to the right to close tube bed.
- ✓ When sensing significant resistance during closing either clamping screws are not open or tube is too large.
- Turn both clamping screws clockwise to clamp tube.
- Check correct mounting by running a short test.

- ✓ The tube may not move during operation.
- To adjust clamping stop rotation, switch off device and open tube bed.
- Turn clamping screws clockwise or counterclockwise to increase or decrease clamping pressure.
- Close tube bed and repeat the test run.

#### **Caution:** **Possible tube damage**

Do not force closing of tube bed. Overtightening leads to damaged tube.

Open clamping screws completely before closing tube bed. Adjust clamping gradually.

Check correct mounting by repeated test runs if necessary.

### ➤ Mount tubes: Multi-channel pump heads

Multi-channel pump heads require cassettes to mount tubes. You will find the right sizes of cassettes in our general catalogue or at the end of this manual.

**Warning:**  
**Risk of crushing!**

While mounting and demounting cassettes you can crush your fingers, if rotation starts accidentally.

Always switch device OFF at main switch before mounting cassettes.

**Caution:**  
**Possible damage of device**

Multi-channel pump heads require a certain combination of cassettes and tubes. Other combinations lead to premature wear or to malfunctions.

Only use combinations of cassettes and tubes fitting your Multi-channel pump head as laid out in this manual (see chapter „Technical data, Table tubes“).

All cassettes operate in both directions. However, we recommend to mount cassettes alternating. This way the clamping can be adjusted more easily.

## Cassettes and tubes: „Cassette small“

„Cassette small“ is suitable for multi-channel pump heads C 4 and C 12 only!


### Caution: Possible damage of device

Stopper on tubes with small inner diameter of 0.8 and 1-7 prevent tubes from moving. Operating tubes without stoppers will lead to malfunction and exclusion from warranty!

Always use stopper tubes when operating these small tubes.  
„Cassette small“ is suitable for two-stopper tubes only!


### Mount tube

- Turn adjustment knob for tube pressure counterclockwise till tube bed is opened to maximum.
- ✓ Wedges are pushed out.
- Insert two-stopper tube into cassette from below.
- Stoppers are placed outside of cassette.
- Pull tube into cassette till stoppers are hold by the side frame.

### Mount cassette

- Click cassette hook onto support rod while holding onto tube.
- Move cassette downwards to click in opposite hook.
- ✓ Tube is positioned between rollers and cassette.
- Push cassette to pump drive.


## Assembly

### Adjust tube pressure

- Turn adjustment knob clockwise till inner side of wedge aligns with mark on adjustment scale between figures 1 and 2.
- ✓ Tube pressure is set.
- Repeat with next cassette.

### Demount cassette

- Push cassette downwards and pull at clamping side to open hook.
- Move cassette upwards while holding onto clamping side to release hook.
- Unhook other side of cassette.

If tube pressure was already adjusted for tube size in previous operations, adjustment of pressure does not have to be repeated.

### Optimize tube pressure

Certain applications require adjustment of tube pressure, for instance to:


- vary flow for certain tube size
- reduce variation in flow caused by system pressure
- extend life time of tubes
- Start the pump.
- Turn adjustment knob counterclockwise to move wedges outwards to higher scale.
- ✓ Tube pressure is reduced to insufficient.
- Slowly adjust by turning knob clockwise till flow and pressure fit your requirements.

### Secure cassettes

Multi-channel pump heads can be operated with reduced number of cassettes. Reduced number of cassettes have to be secured to prevent moving.

Both support rods are equipped with an O-ring to prevent cassettes from moving.

- Push O-rings on both sides away from pump head.
- Mount required number of cassettes.
- Push O-Ring on both sides against the front cassette.
- ✓ Cassettes are secured.


## Tube extension

Optional PTFE tube couplings are available for extension. Bulk hose goods can be connected with these couplings.

### Components:

- 2x Sleeve nuts
- 1x Screw joint
- 2x3 Pcs. O-rings
- 2x Ferrules

### Tube sizes:


#### 0.2; 0.5; 0.9 and 1.4

- Insert tubes into sleeve nuts.
- Push 3 O-rings over each tube.
- Screw sleeve nuts hand-tight onto screw joints.

#### Tube size 2.8

- Trim tube ends diagonally.
- Insert tubes into sleeve nuts.
- Trim tube ends horizontally.
- Push ferrules with flat side onto tube.
- Insert tubes with ferrule tip into screw joints at both sides.
- Screw sleeve nuts hand-tight onto screw joints.

## Cassettes and tubes: „Cassette medium“

„Cassette medium“ is suitable for multi-channel pump head C 8 only!

- Mount suitable tubes to cassettes. (see chapter „Technical data, Table tubes“).
- Mount Cassettes as described under “Cassette small”.
- Pull slightly at both ends of tube to prevent an eventual loop.
- Set tube pressure via adjustment knob between figures 3 and 4 of scale.



# Cassettes and tubes: „Cassette large“

„Cassette large“ is suitable for multi-channel pump head C 8 only!

- Mount suitable tubes to cassettes. (see chapter „Technical data, Table tubes“).
- Set tube size by moving slider to fitting size.
- Mount Cassettes as described under „Cassette small“.
- Pull slightly at both ends of tube to prevent an eventual loop.
- Press slider slightly inwards.
- Position the slider marking at the prominent white lines according to your tube size.
  - Tube size 1.7, 3.1 to marking 16/14
  - Tube size 4.8 to marking 25
  - Tube size 6.3 to marking 17
- Turn adjustment knob clockwise till inner side of wedge aligns with mark on adjustment scale between figures 3 and 4.

### ➤ Mount several pump

When operating several pumps simultaneously you can place up to 2 pumps on top of each other.


#### **Warning:** **Risk of slipping device**

In operation the device on top may start moving and slide off the set up area.

The set up area must be clean, be properly leveled and the devices must be positioned exactly on top of each other.

- Position pump on top of each other in a way that left and right sides as well as back sides are aligned.

### ➤ Dismantling, transportation and storage

#### Dismantling

#### **Warning:** **Risk of injury!**

While the device is plugged in you might be at harm if the device is accidentally switched to ON mode.

Prior to dismantling device ensure that the device is switched off and the plug is pulled from the electric socket.

- ✓ The device is switched off.
- Disconnect plug from power socket.
- Remove all vessels.
- Disassemble attachments.

#### Transportation and storage

- Store and transport the device and its components only if they were emptied, cleaned and dismantled.
- To carry the device wear safety gloves and lift the unit from below.
- Store and transport the device and parts of device in the original packing material or alternatively in an adequate container to prevent damage.
- Seal the package carefully against unauthorized or accidental opening.
- Store the device in a dry and frost-free place.

Improper storage and transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

## ➤ Scope of delivery

Component	Variant	Quantity	Product number
Hei-FLOW	Value 01*	1	523-50010-00
	Value 06*	1	523-50060-00
	Advantage 01*	1	523-51010-00
	Advantage 06*	1	523-51060-00
	Precision 01*	1	523-52010-00
	Precision 06*	1	523-52060-00
Operating Manual English / German		1	01-005-002-75
Warranty registration / Confirmation of condition		1	01-006-002-78
EC Declaration of Conformity		1	01-001-025-06

\* The product numbers are valid for 230 V devices in the european community. Please contact Heidolph Instruments for further product numbers on country specific varieties.

## ➤ Accessories

Fitting variant	Component	WT*	Quantity	Product number
all pump heads	Pump head SP standard	1.6 mm	1	523-43010-00
	Pump head SP standard	2.5 mm	1	523-43030-00
	Pump head SP vario	1.6 mm - 2.5 mm	1	523-45110-00
	Pump head SP quick	1.6 mm	1	527-11100-00
	Pump head SP quick	2.5 mm	1	527-11300-00
Only Hei-FLOW: ▪ Value 01 ▪ Advantage 01 ▪ Precision 01	Pump head for 4 x Casette small C 4	0.9 mm	1	524-80420-00
	Pump head for 8 x Casette medium C 8	1.6 mm	1	524-80810-00
		for 4 x Casette large	1.6 mm	
	Pump head for 12 x Casette small C 12	0.9 mm	1	524-81220-00

\* WT = wall thickness

Additional accessories are shown in the general catalog or on our website at [www.heidolph-instruments.com](http://www.heidolph-instruments.com).

## ➤ Technical data

### Hei-FLOW Peristaltic pumps

Standard supply voltage (V)	110-240 V, 50/60 Hz
Power consumption (W)	100
Protective class (IEC 61140)	1 ⊕
Protection class (IEC 60529)	IP 55
Sound pressure level (dB(A)) (based on IEC 61010)	< 55
Motor	DC-motor
Motor power (W)	100
Speed control	digital
Speed accuracy %	± 0,5
Operating mode	continuous 24 hours / 7 days
Rotating direction	right / left
Dimensions (excl. pump head) (w x h x d) (mm)	166 x 225 x 256

Model	Value		Advantage		Precision	
	01	06	01	06	01	06
Overheat protection	yes		yes, electronic		current limiter	
Display speed	analog		analog		digital	
Speed range	10-120	50-600	5-120	24-600	5-120	24-600
Flow rate single-channel pump head (ml/min)	0.85-861	4.0-4151	0.38-813	2.0-4056	0,36-813	2.0-4056
Flow rate multi-channel pump head (ml/min)	0.005-364	-	0.005-329	-	0.005-329	-
Accuracy flow*(%)	± 5		± 3,5		± 1 ± 2	
Display flow	-		-		digital	
Volume dosing (ml)	-		-		0.001-9999	
Intermittant dosing (ml)	-		-		0.001-9999 pauses 0.1 sec - 750 h	
Smooth start	-		-		yes	
Elektronic break	-		-		yes	
Interface digital	-		-		RS 232	
Interface analog	-		(0-10 V / 4-20 mA) speed, direction, start/stop			
Connection for foot pedal	-		yes			
Weight (excl. pump head) (kg)	7.6	7.1	7.6	7.3	7.7	7.3

\* Flow rates pertain to water, no counter pressure

## Table tubes

Flow rates are guidelines and pertain to Tygon® (standard) tube, media water.

Single-channel pumps		Wall thickness 1,6 mm				
Ø inner (mm)		0.8	1.7	3.1	4.8	6.4
Ø outer (mm)		4	4.9	6.3	8	9.5
Pressure max. continuous (bar)		0.7			0.5	
Pressure max. short time (bar)		1.7			1.5	
Suction height (mWs)		8.8			8.8	6.7
Pump head	Pump	Flow rate H <sub>2</sub> O (ml/min)				
SP standard SP vario (Pu 2)	Hei-FLOW Value 01	-	5-61	19-223	44-519	75-861
	Hei-FLOW Advantage 01, Hei-FLOW Precision 01	-	2-55	9-221	21-530	33-813
	Hei-FLOW Value 06	-	22-249	93-1,037	228-2,613	364-4,151
	Hei-FLOW Advantage 06, Hei-FLOW Precision 06	-	11-257	43-1,017	105-2,549	167-4,056
SP quick (Pu 1)	Hei-FLOW Value 01	0.83-9	3-41	11-134	25-292	36-413
	Hei-FLOW Advantage 01, Hei-FLOW Precision 01	0.38-9	2-40	5-126	12-233	17-409
	Hei-FLOW Value 06	4-35	17-197	57-695	123-1,494	186-1,765
	Hei-FLOW Advantage 06, Hei-FLOW Precision 06	2-33	8-186	26-653	59-1,529	89-2,072

Single-channel pumps		Wall thickness 2,5 mm		
Ø inner (mm)		4.8	6.4	7.9
Ø outer (mm)		9.8	11.3	12.9
Pressure max. continuous (bar)		0.8		
Pressure max. short time (bar)		1.8		
Suction height (mWs)		8.8		
Pump head	Pump	Flow rate H <sub>2</sub> O (ml/min)		
SP standard SP vario (Pu 2)	Hei-FLOW Value 01	42-493	68-773	-
	Hei-FLOW Advantage 01, Hei-FLOW Precision 01	15-491	28-769	-
	Hei-FLOW Value 06	203-2,426	313-3,782	-
	Hei-FLOW Advantage 06, Hei-FLOW Precision 06	92-2,390	139-3,821	-
SP quick (Pu 1)	Hei-FLOW Value 01	26-299	38-454	50-636
	Hei-FLOW Advantage 01, Hei-FLOW Precision 01	12-299	18-435	25-630
	Hei-FLOW Value 06	123-1,580	180-2,411	257-3,436
	Hei-FLOW Advantage 06, Hei-FLOW Precision 06	58-1,527	85-2,248	113-3,171

## Attachments

### Multi-channel pumps

### Wall thickness 0.9

Ø inner (mm)		0.25	0.51	0.89	1.42	2.79	
Ø outer (mm)		2.05	2.31	2.69	3.22	4.59	
Pressure max. continuous (bar)				0.5			
Pressure max. short time (bar)				1.5			
Suction height (mWs)				7			
Pump	Pump head	Cassette	Flow rate H <sub>2</sub> O (ml/min)				
Hei-FLOW Value 01	C 4: (Pu 9)	small: max. 4 Cassettes	0.04-0.53	0.17-2	0.57-6	1-15	4-37
	C 12: (Pu 8)	small: max. 12 Cassettes	0.005-0.11	0.02-0.42	0.10-1	0.23-3	0.69-8
Hei-FLOW Advantage 01, Hei-FLOW Precision 01	C 4: (Pu 9)	small: max. 4 Cassettes	0.02-0.49	0.08-2	0.24-6	0.60-14	2-36
	C 12: (Pu 8)	small: max. 12 Cassettes	0.005-0.11	0.01-0.54	0.03-1	0.10-3	0.29-9

### Multi-channel pumps

### Wall thickness 1.6

Ø inner (mm)		0.8	1.7	3.1	4.8	6.4	
Ø outer (mm)		4	4.9	6.3	8	9.5	
Pressure max. continuous (bar)		0.7		0.7		0.5	
Pressure max. short time (bar)		1.7		1.7		1.5	
Suction height (mWs)		7		7		6.7	
Pump	Pump head	Cassette	Flow rate H <sub>2</sub> O (ml/min)				
Hei-FLOW Value 01	C 8: (Pu 0)	medium: max. 8 Cassettes	0.55-6.97	2.17-27	-	-	-
		large: max. 4 Cassettes	-	2-27	7-85	18-246	26-364
Hei-FLOW Advantage 01, Hei-FLOW Precision 01	C 8: (Pu 0)	medium: max. 8 Cassettes	0.24-7	1-26	-	-	-
		large: max. 4 Cassettes	-	1-27	4-90	8-192	11-329

For high dosing accuracy choose a tube with small inner diameter.  
For fast dispensing choose a tube with large inner diameter.


Flow rate depends on tube material, dispensed media and temperature. Flow rate can therefore vary and may deviate significantly from indicated values.


# Performance range pump heads

## Single-channel pumps:

### SP standard, SP vario


### SP quick


## Multi-channel pumps:


### C 4

### C 12

### C 8

## ➤ Contact / Technical Service

### Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

### **Warning!** **Danger of poisoning!**

Contaminated devices can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
- take adequate protective action towards our parts receiving and service personnel
- mark the pack in accordance with Ordinance on Hazardous Substances

A "Confirmation of Condition" form can be found at the end of this operating manual.

→ Prior to shipping the device for repair, complete a copy of this form and submit it in advance.

### Contact details:

#### **Heidolph Instruments Germany**

Heidolph Instruments GmbH & Co. KG  
Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach / Germany

Tel.: +49 – 9122 - 9920-74  
Fax: +49 – 9122 - 9920-84

E-Mail: [service@heidolph.de](mailto:service@heidolph.de)  
[www.heidolph.com](http://www.heidolph.com)

#### **Heidolph Instruments North America**

Phone: 1-866-650-9604  
E-mail: [service@heidolph.com](mailto:service@heidolph.com)  
[www.heidolphNA.com](http://www.heidolphNA.com)

#### **Heidolph Instruments United Kingdom**

Phone: 01799 - 5133-20  
E-mail: [service@radleys.co.uk](mailto:service@radleys.co.uk)  
[www.heidolph-instruments.co.uk](http://www.heidolph-instruments.co.uk)

#### **All other countries**

You will find contact details of your local Heidolph distributor at  
[www.heidolph.com](http://www.heidolph.com)

## ➤ Warranty

---

Heidolph Instruments provides a 3 year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online ([www.heidolph.com](http://www.heidolph.com)). The warranty period begins with the date of registration. The serial number is also valid without registration.

This warranty covers material and manufacturing defects.

---

In the event of a material or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damage incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please notify Heidolph Instruments or your local Heidolph distributor if you wish to make a warranty claim.

## ➤ Confirmation of condition

→ In the case of repair, copy and complete the Confirmation of condition form and send it to Heidolph Instruments.

### 1. Details about the device

Product number \_\_\_\_\_

Serial number \_\_\_\_\_

Reason for repair \_\_\_\_\_

\_\_\_\_\_

### 2. Has the device been cleaned, decontaminated/sterilized?

Yes \_\_\_\_\_ No \_\_\_\_\_

### 3. Is the device in a condition which does not represent any health threats for the staff of our service department?

Yes \_\_\_\_\_ No \_\_\_\_\_

If not, which substances has the device come into contact with?  
 \_\_\_\_\_  
 \_\_\_\_\_

### 4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

<b>Date</b>	<b>Signature</b>
-------------	------------------

### Company stamp

### Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG