



# Cultivation Systems From Discovery to Production




Engineering for Life

# micro-Matrix |

## 24 Bioreactors in a Convenient Microtiter Format

- Do more cultures in less time
- Liquid addition to the 24 wells
- PID controls of individual wells
- User friendly software enables more configuration for more complex recipes
- Small footprint on the bench in the lab




The unique micro-Matrix offers total control over 24 independent bioreactors in a simple microtiter plate footprint. Each of the 24 bioreactors on a plate offers independent controls, like its larger stirred-tank relatives:

- pH (measurement and two-sided control)
- temperature (measurement and two-sided control, including plate-wide gradients)
- dissolved oxygen control (measurement and two-sided control)
- individual liquid additions (including feeding profiles)
- up to 4 separate gas additions per row (individually controlled)

The micro-Matrix is a true scale-down of small scale bioreactors. The bioreactor square well cassette design is based upon the SBS-format microtiter plates that maximize mixing, optimize gas transfer, and seamlessly integrate into lab automation protocols. The PC-based human interface of the micro-Matrix reflects our popular my-Control interface and offers simple, intuitive interaction with each of the 24 bioreactors. Integrated LEDs indicate the status of the bioreactors (inactive / active / alarm) with color-based feedback, so that operators can get instant process information at one quick glance.

## Features

- Individual monitor and control: pH, dissolved oxygen, temperature
- Liquid feeding enables true scale-down studies
- Simple, powerful software streamlines operator workflow
- Operating 24 bioreactors in parallel
- Screening tool



## Applications

- Screening of cell-line and microbial cultures
- Process development studies
- Process optimization studies
- Small volume cultivations
- Anaerobic conditions for small volumes



## Software

The micro-Matrix software offers an easy way to operate 24 bioreactors in parallel, plus simple comparisons of large numbers of experimental cultures. Using an instinctive left-to-right progression, the interface guides operators through instrument configuration, control strategy definition, experimental setup, and

data visualization. It is also possible to define time- and event-based control actions, and all data can be exported from the instrument mid-cultivation. At any moment in the run, an Excel file can be generated with graphs of all measurements.



## Measurement and control

Each bioreactor has its own PID controller for pH, dissolved oxygen, and temperature. Individual pH control can be achieved via gas addition, liquid addition, or a combination of both. Dissolved oxygen level can be individually controlled by up to four gas additions per bioreactor. Temperature can be individually controlled by the integrated cooling and heating system on a


well-by-well basis, and users can define temperature gradients across the set of bioreactors. Advanced control strategies are also available to users interested in cascade controls, time-based setpoint changes, and event-triggered liquid feeds. Liquid feeding is also available using varying types of addition profiles.

## Technical data

Each of the four gas feeds and one liquid entry are controlled individually per bioreactor (5 additions per reactor).

- Independent gas supplies: for Air / O<sub>2</sub> / N<sub>2</sub> / CO<sub>2</sub> / NH<sub>3</sub>\*
- Integrated liquid feed individually per well
- Cascade up to five actuators per control loop per bioreactor
- Control strategy is user definable per bioreactor per experiment
- Temperature measurement and control per reactor
- pH measurement and control per reactor
- DO measurement and control per reactor
- Orbital shaker 0-400 rpm @25 mm orbit

\*Ammonia gas is supplied via Applikon designed ammonia pressure vessel for safety conditions.



Schematic diagram illustrating the functioning of the micro-Matrix control loops.

# Specifications

<b>Total volume per bioreactor well</b>	10 mL	
<b>Working volume per bioreactor well</b>	2 - 5 mL	
<b>Orbit diameter</b>	25 mm	
<b>Rpm range</b>	0 - 400 rpm	
<b>Gas delivery</b>	0.1 - 15 mL / min	Per gas, per well
Number of gasses	4	
Gas feed pressure	> 1.0 (extern)	
<b>Liquid delivery</b>	200 nL	per bolus (approx.)
Liquid feed pressure	> 2.0 bar (extern)	
<b>Temperature measurement range</b>	0 - 45 °C	
Resolution	0.1 °C	
Accuracy	0.5 °C	
Temperature control range	15 - 45 °C	
ΔT between adjacent wells	1 °C	
Heat-up time	1 °C / min	20 °C to 37 °C
Control accuracy	± 0.1 °C	
<b>pH measurement</b>	5.5 - 8.5 pH	
Resolution	< 0.05 pH	@ pH 7.2
Accuracy	< 0.1 pH	@ pH 7.2
pH control range	5.8 - 8.2 pH	
Control accuracy	± 0.1 pH	
<b>DO measurement</b>	0 - 150 %	Air saturation
Resolution	< 1 %	Air saturation
Accuracy	< 0.5%	@ 0% air saturation
	< 3%	@ 100% air saturation
DO control range	0 - 150 %	
Control accuracy	± 5%	Air saturation
<b>Dimensions</b>	415 mm	Height
	570 mm	Depth
	560 mm	Width
<b>Weight</b>	65 kg	
<b>Oxygen Transfer Rate</b>	250 mmol / l / h	
<b>Mixing times</b>	< 1 sec	


## Connections

The micro-Matrix requires connection to a standard 230/115V AC power supply and pressurized gas and (air, oxygen, nitrogen, carbon dioxide). The micro-Matrix is connected to a computer through a standard Ethernet port.

## Related products

- The polystyrene bioreactor cassettes conform to the SBS standard for 24 well culture plates (128 x 86 mm). The total volume per bioreactor is 10 mL with optimal working volumes from 2 to 5 mL. The square bioreactors are designed based upon our popular 24-well microtiter plate design to optimize mixing and gas transfer, and Applikon has characterized these plates to offer advice on proper selection of working volume and agitation to achieve your process goals. Cassettes are delivered beta-sterilized and sealed in a light shielding package.
- Filter bar





## America

### Argentina

AP Biotech  
Phone: +54 11 4292 1611  
[www.apbiotech.com.ar](http://www.apbiotech.com.ar)

### Brasil

Lobov Cientifica  
Phone: +55 11 3829-8040  
[www.lobov.com.br](http://www.lobov.com.br)

### Mexico

Applikon Biotechnology Mexico  
Phone: +52 55 5868 7814  
[www.applikon-biotechnology.com](http://www.applikon-biotechnology.com)

### Peru

Applikon Biotechnology Peru  
Phone: +51 1657 1230  
[www.applikon-biotechnology.com](http://www.applikon-biotechnology.com)

## Europe

### Belgium

Applikon Biotechnology BV  
Phone: +31 10 208 3555  
[www.applikon-biotechnology.com](http://www.applikon-biotechnology.com)

### Czech Republic

Vekamaf spol. s.r.o.  
Phone: +421 907 753 338  
[www.vekamaf.cz](http://www.vekamaf.cz)

### Denmark

Holm & Halby A/S  
Phone: +45 432 694 00  
[www.holm-halby.dk](http://www.holm-halby.dk)

### Estonia

Biotecha Eesti OÜ  
Tel. 00372 659 7101

### France

System-c bioprocess  
Phone: +33 475 548 600  
[www.system-c-bioprocess.com](http://www.system-c-bioprocess.com)

### Germany & Austria

I&L Biosystems GmbH  
Phone: +49 222 391 920  
[www.il-biosystems.com](http://www.il-biosystems.com)

### Hungary

Holimex Mérnöki és Kereskedelmi Kft.  
Phone: +36 1 391 44 00  
[www.holimex.hu](http://www.holimex.hu)

### Ireland

Asistec  
RTD Technology Ltd,  
Phone: +353 1 281 4650 / 281 4064  
[www.asistec.ie](http://www.asistec.ie)

### Italy

Steroglass S.r.l.  
Phone: +39 075 609 091  
[www.steroglass.it](http://www.steroglass.it)

### Lithuania

BIOTECHA UAB  
Phone: +370 5 237 6007  
[www.biotecha.lt](http://www.biotecha.lt)

### Poland

Labo Baza  
Phone: +48 61 812 57 45  
[www.labobaza.pl](http://www.labobaza.pl)

### Russia

ABTEK Bioprocess Solutions  
Phone: +7 (495) 723 6497  
[www.abtek.ru](http://www.abtek.ru)

### Slovak Republic

Vekamaf spol. s.r.o.  
Phone: +421 907 753 338  
[www.vekamaf.cz](http://www.vekamaf.cz)

### Spain & Portugal

VERTEX Technics, S.L.  
Phone: +34 93 223 33 33  
[www.vertex.es](http://www.vertex.es)

### Sweden

AB Ninolab  
Phone: +46 8 590 860 70  
[www.ninolab.se](http://www.ninolab.se)

### Switzerland

ReseaChem GmbH  
Phone: +41 344 240 310  
[www.reseachem.ch](http://www.reseachem.ch)

### Turkey

Ant Teknik Cihazlar Ltd. Sti.  
Phone: +90 216 422 67 00  
[www.antteknik.com](http://www.antteknik.com)

## Australia & New Zealand

### Australia

John Morris Scientific  
Phone: +61 2 9496 4200  
[www.johnmorris.com.au](http://www.johnmorris.com.au)

## Asia

### China

Applitech Pharma  
Phone: +86 203 482 1111 (ext. 8216)  
[www.applitechpharma.com](http://www.applitechpharma.com)

### India

Spinco Biotech Pvt. Ltd.  
Phone: +91 44 2434 0174  
[www.spincootech.com](http://www.spincootech.com)

### Indonesia

Pt. Arico Sainsindo Nusantara  
Phone: +62 21 580 6342 / 3  
[www.ptarico.com](http://www.ptarico.com)

### Japan

Sanyo Trading Co., Ltd.  
Phone: +81 3 3518 1111  
[www.sanyo-trading.co.jp](http://www.sanyo-trading.co.jp)

### Korea

Seoulin Bioscience Co., Ltd.  
Phone: +82 1670 5911  
[www.seoulin.co.kr](http://www.seoulin.co.kr)

### Singapore

Zil Sci Pte. Ltd.  
Phone: +65 676 729 29  
[www.zilsci.com](http://www.zilsci.com)

### Taiwan

Tseng Hsiang Life Science Ltd.  
Phone: +86 2 2785 1156  
[www.thco.com.tw](http://www.thco.com.tw)

### Thailand

Meditop Co., Ltd  
Phone: +66 2933 1133  
[www.meditopthailand.com](http://www.meditopthailand.com)

## Middle East & Africa

### Israel

Medent LTD  
Phone: +972 8 623 7984  
[www.medent.co.il](http://www.medent.co.il)

### South Africa

Separations  
Phone: +27 11 919 1000  
[www.separations.co.za](http://www.separations.co.za)

### Headquarters The Netherlands

**Applikon Biotechnology BV**  
Heertjeslaan 2  
2629 JG Delft  
The Netherlands  
T. +31(0)10-208 35 55  
F. +31(0)10-208 35 05  
E. [info@applikon-biotechnology.com](mailto:info@applikon-biotechnology.com)

### USA Main Westcoast Office

**Applikon Biotechnology Inc.**  
1180 Chess Drive  
Foster City CA 94404  
Phone: +1 650 578 1396

### Latin America Office

Phone: +52 55 5868 7814

### United Kingdom Office

**Applikon Biotechnology**  
Basepoint Business Centre  
Oakfield Close  
Tewkesbury Business Park,  
Tewkesbury, Gloucestershire, GL20 8SD  
Phone: +44 1684 851 281  
United Kingdom



Engineering for Life  
[www.applikon-biotechnology.com](http://www.applikon-biotechnology.com)