



One

SINGLE WALL
FERMENTERS/
BIOREACTORS



SOLARIS
BIOTECH SOLUTIONS

SINGLE WALL FERMENTERS/BIOREACTORS

one

The **ONE** platform is ultimate Solaris entry level system. It offers multiple autoclavable vessel sizes and designs from 2 up to 10 L total volume. Various aspect ratios and thermoregulation designs are also available. The system is 100% configurable, built with high quality components, and offered at a competitive price.

ONE typical applications includes the following:

- Education & Basic research
- Scale-up and scale-down studies
- Process development and optimization

One can be used for:

- Biopharmaceutical
- Biofuels
- Food industry
- Bioremediation
- Bioplastic
- Cosmeceutical
- Nutraceutical



**WHY TO
INVEST
IN THIS PRODUCT**

The essential
solution
shaped on
your needs

Entry Level
Fully
configurable


SINGLE WALL FERMENTERS/BIOREACTORS

one

Benefits

LEONARDO software license: smart controller designed to provide an high level of automated management of the fermentation/ cultivation processes
Batch, Fed batch or continuous processes

Different gas mixing strategies with up to 5 TMFC or up to 4 rotometers



After sales assistance

Feeding options: single or multi feed connector



Powerful/ Accurate **brushless motor**, from 1 to 2000 RPM.
Online absorbed Torques (Nm) and Power (W) measurements obtaining an indirect density indication of the culture broth.

pH, dO2, antifoam, level and temperature controls available



LEDA safe sterile sampling system

Safety: pressure relief valve included in each unit

Stainless-steel cube PCS



Up to n. 4 WM 114 Watson Marlow fixed speed pumps available

5 different volumes and 2 different ratio H/D

Single wall vessel, with heating blanket and cooling finger (available options)



Modbus Digital sensors

Why a digital sensor?

Digital sensors (pH, dO2) have been integrated to the Solaris PCS and Leonardo controlling software, giving the user many benefits over traditional analog sensor outputs. Such benefits include a robust communication protocol not susceptible to signal loss, in-software sensor diagnostic information, parallel calibration/batch calibrations and more.



Sensor life traceability

Reducing background noise

GAS MIXING

Configure your system with a selectable number of TMFCs (up to 5) and/or rotameters (up to 4).

- Various agitator and baffle designs available
- Automatic gas mixing algorithms
- Toro, sintered and other spargers available



LEDA sterile sampling system

Technical specifications

Material	VALOX resin (external) silicone (internal)
Autoclavable	121-133°C (up to 30 minutes)
Residual volume	0.04 mL
Flow rate	165 mL/minute



- Sterile single use sampling system up to 180 sterile sampling per batch.
- Needlefree connector is designed to reduce the risk of contamination during sampling.
- The sterile combination of a syringe (3-5-10-30 ml) and a non return valve guarantees the sterility after sampling until the next use

Smart PCS

- Compact PCS
- 35 x 35 x 35 cm
- Lean design
- Software license on USB flash key



Leonardo 3.0

USER-FRIENDLY SOFTWARE

Solaris controlling software offers a simply laid out, yet powerful platform for experimental design planning and process control. The graphical user interface enables the intuitive selection and adjustment of control functions. Extracted data is compatible with Window Excel but, in addition, Solaris offers a platform where fermentation data can be easily exported in real time and thus managed. This software license is included in the supply and can be installed on an unlimited number of the client's PC or laptops.



Features:

- Home with Multi-level password protection
- Workflow settings (before cultivation, cultivation, custom phases)
- Synoptic page with manual operation of all the actuators (pumps, valves etc.)
- Continuous trend graphs representation to track, print and export data on up to 6 processes and set point variables.
- Different dynamic zooms and configurations in a time frame that can be set interactively
- Cascade control (optional) and profile programs
- Pumps configuration and calibration
- On line parameters calibration
- PID setting
- After sale assistance
- Possibility of saving up to recipes for repeat usage
- Print-out of hard copy of each screen

Data sheet

Vessel					
Solaris Code	One 2.0	One 4.0	One 6.5	One 8.0	One 10.0
Production Code	onest2.0	onest4.0	onest6.5	onest8.0	onest10.0
Total Volume (L)	2,00	4,00	6,50	8,00	10,00
Ratio D/H	1:3,0	1:3,0	1:2,5	1:3,0	1:3,0
Min. Working Volume (L)	0,35	0,60	1,10	1,10	1,60
Max. Working Volume (L)	1,40	2,80	4,50	5,50	7,0
Max. temperature	70°C				
Operating pressure	< 0.5 bar				
Headplate Ports (n.10 in One 2.0; n.13 in the others)	10: n.1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Temperature, n.1 Multifeed. n.2 Sensors DN12, n.1 Spare. 13: n.1 Agitation Group, n.1 Gas Sparger, n.1 Gas Overlay, n.1 Gas Out/Condenser, n.1 Sampling/Harvesting, n.1 Sterile Sampling System, n.1 Temperature, n.1 Multifeed, n.2 Sensors DN12, n.3 Spare.				
Design	Borosilicate glass vessel (single wall)				
Materials	Borosilicate Glass and AISI 316 L				
Sensors length (mm)					
pH	325	425	425	425	425
dO ₂	325	425	425	425	425
Dimensions for autoclave (with Condenser)					
Height (mm)	610	705	705	790	790
Diameter (mm)	275	285	315	315	335
Stirring					
Drive	Brushless Motor				
Speed (rpm)	1-1900	1-1800	1-1700	1-1700	1-1700
Nominal Torque (Nm)	0,9	0,9	0,9	1,1	1,1
Impellers	Select from: Rushtons impellers, Marine Impellers, Pitched blade				
Thermoregulation					
Control	PID Control - Accuracy 0,1 °C - Jacketed with n. 2 Electric Cartridge Heaters and cooling valve				
Total Heater Power (W)	400	600	700	700	700
Gas Control & Gas Mixing					
Sparger and overlay Gas Control	TMFC				
Gas Mixing (Air,CO ₂ ,O ₂ ,N ₂)	n.1 TMFC (included in entry level) + n.4 solenoid valves or + n. of additional TMFC (up to n.4)				
Sparger type	Select from: Toro type (ring), sintered microbubbling - both provided with 0,22 µm sintered filter				
Gas Out	n. 1 Condenser + 0,22 µm sinterized filter				
Peristaltic Pumps					
	up to n.4 Watson Marlow type 114, fixed speed, max. 60 rpm, volumetric flow 0,5-51 ml/min, function assignable from software				
Controller					
Master Control Module	35 x 35 x 35 cm				
HMI with Leonardo software	Licence				

Controls

Temperature	
Sensor	PT100
Accuracy	0,1 °C
Control system	Measuring resident in Leonardo 3.0 software
Control range	0 - 70°C
pH	
Sensor	Digital sensor
Sensitivity	57 to 59 mV/pH
Control system	Measuring resident in Leonardo 3.0 software
Control range	0 - 14
Operation temperature	0 - 130°C
Pressure range	0 - 6 bar
dO ₂	
Sensor	Digital Optical sensor
Accuracy	±0.05%-vol, 21±0.2%-vol, 50±0.5%-vol
Control system	Measuring resident in Leonardo 3.0 software
Control range	0,05 - 300% air saturation
Operation temperature	-10 - 130°C
Pressure range	0 - 12 bar
Antifoam/Level	
Sensor	Solaris sensor
Control	Measuring resident in Leonardo 3.0 software

INTEGRATED IN THE PCS


Chiller

- Optionally ONE can be equipped with a chiller for heat removal from your culture minimizing lab water usage
- Using this system you don't need a water supply line in your lab
- Cost-effective cooling of fermenters
- Easy operation
- Refregerant level monitoring



Chiller data sheet

Working temperature range	-10°C / +40°C
Temperature stability	±0.5
Power consumption	0.7 kW
Filling volume range	2-8 L
Cooling output at 20°C measured with ethanol	0.25-0.60 kW
Cooling output at 10°C measured with ethanol	0.20-0.50 kW
Cooling output at 0°C measured with ethanol	0.15-0.36 kW
Cooling output at -10°C measured with ethanol	0.09-0.15 kW
Pump pressure max.	0.35-1.30 bar
Pump flow max.	16-35 L/min.



SOLARIS

BIOTECH SOLUTIONS

SOLARIS BIOTECHNOLOGY srl
Via Bachelet, 58 - 46047 Porto Mantovano
Mantova - Italy
Phone: +39 0376 408760
Fax: +39 0376 385108
Email: info@solarisbiotech.com
www.solarisbiotech.com