

Mini-Pipette washer MPS



Washing procedure

1. Place the contamined pipettes into

2. Choose the washing program and


the rinser by using the basket system

Product description

The MPS is an fully automatic pipette washer with two different washing programs for all common kinds of glass pipettes. The MPS has one water connection and can be used with demineralized water or with normal tap water.

If the MPS will be used with demineralized water the dosage of liquid detergent takes place automatically. When using tap water the filling in of detergent has to be done manually. A buffer solution will be added automatically to the water tap.

- \oplus small dimensions for little laboratories
- \oplus two washing programs
- \oplus operation with demineralized <u>or</u> tap water
- ⊕ program runs automatically
- suitable for all kinds of common glass pipettes
- ⊕ high cleaning capacity
- \oplus surface considerate cleaning process



Parameters for:

3. Fully automatic program run

press Start button 1 or 2

4. The clean and dry pipettes can be Taken out of the machine

GEWO Feinmechanik GmbH

Technical features

Connection requirements

mains:	230V / 10 A
water:	demineralised <u>or</u> tap water
	0,5 to 6 bar water pressure
	max. 60° C water inlet temperature

Consumption

powert:4 kWhwater:max. 55 ldetergent:0 to 40 ml

Temperature

cooking: 70°C or 95°C drying: 130°C

Programs for on kind of water

number:	2
duration:	Program 1 (Standard): 2 h 50 min
	Program 2 (Intensive): 3 h 50 min

Capacity: mind. 1 basket Ø 135 mm max. 4 basket Ø 60,3 mm

Weight/ material: 40 kg / stainless steel housing

Dimensions:

breadth:	250 mm
depth:	485 mm
height:	920 mm (device top edge with closed cover plate)
height:	1150 mm (with open cover plate)



H.H

Labortechnik für Forschung und Industrie

Abteilung

ZEL

Included accessories

1 pcs. basket Ø135 x 300 mm

1 liter liquid detergent

Technical documentation and connection material



Pipette capacity

Measurement pipettes									
Volume	1,0 ml	2,0 ml	5,0 ml	10,0 ml	20,0 ml	50,0 ml	100,0 ml		
Rinser	180	160	145	75	55	25	12		
Full pipettes									
Full pipett	es								
Full pipett Volume	es 2,0 ml	3,0 ml	5,0 ml	10,0 ml	20,0 ml	50,0 ml	100,0 ml		