



# SMART 6

MOISTURE/SOLIDS ANALYZER  
One System, Limitless Testing



### Use At-Line

An IP65 compliant option that provides protection from water and the elements

### Intelli-Temp™

A patented technique that offers constant monitoring of weight and power regulation, while providing optimum repeatability

### Statistical Analysis

Easily review data with options to perform statistical analyses, directly through the system

### Octawave™ Cavity

A patented design created for uniform heating

### iPower®

The patented technique that allows samples of any moisture content level to be dried rapidly

### Interactive Touchscreen

High resolution touchscreen interface with menu-driven software  
Intuitive and simple test process that clearly displays results and includes training videos

### Data Output

0.01% resolution  
0.01% to 99.99% range

### Balance

25 g capacity  
0.1 mg resolution

### Data Storage

Unlimited data storage  
Easy data export

### Exhaust

Safely remove volatiles  
free up hood space

### Accessory Ports

4 x USB-a, 2 x USB-b, 1 x ethernet

### Printer Options

Optional internal thermal printer  
(request list of compliant external printers)

### Electrical requirements

115 VAC  $\pm 10\%$ , 60 Hz, 15 Amps  
220-240 VAC  $\pm 10\%$ , 50 Hz, 10 Amps  
100 - 200 VAC  $\pm 10\%$ , 50/60 Hz, 15 Amps




## SMART 6 Moisture/Solids Analyzer

One System, Limitless Testing.

Once again, CEM has revolutionized moisture/volatile analysis for every application and industry with the SMART 6™. The SMART 6 features two major innovations, compared to other moisture analyzers: the ability to accurately analyze both wet and dry samples, and a significant reduction in sample drying times. These benefits are possible through the development of iPower® technology, which features a patented multi-frequency energy source for sample drying. Drying time now up to 10x faster than IR drying systems, which have traditionally served some markets. Never before has a moisture analyzer been able to offer so much, with the ability to have **the fastest testing and most accurate data for any sample type.**





## Test Any Sample, Faster - The iPower Advantage

iPower is an innovative development in moisture analysis for the SMART 6 that utilizes dual electromagnetic wavelengths, controlled by an intelligent processing system. iPower is able to apply the correct level of energy to the sample at the optimum time, using low frequencies to quickly evaporate free moisture and higher frequencies to carefully remove any water molecules bound to the sample structure, resulting in **unparalleled precision and speed**.

|                        | iPower | NIR | IR | MW |
|------------------------|--------|-----|----|----|
| Direct Analysis        | ✓      |     | ✓  | ✓  |
| Rapid Analysis         | ✓      | ✓   |    |    |
| Low Moisture           | ✓      | ✓   | ✓  |    |
| High Moisture          | ✓      | ✓   |    | ✓  |
| No Costly Calibrations | ✓      |     | ✓  | ✓  |





## Speed

up to 10x faster than IR and Halogen analyzers

iDri™ analysis for results in as little as 90 seconds

Library of methods available for less time spent on method development and validation



## Accuracy

AOAC approved technology

Certified as more accurate than IR, Halogen, or NIR by numerous accredited and unbiased third parties

New method controls improve accuracy across a range of products when coupled with iPower



## Versatility

iPower now gives the SMART 6 the ability to test any product, wet or dry, polar or nonpolar

Now only 1 analyzer is needed for any lab

Multiple data management options fit your exact laboratory needs

|                      | SMART 6    |         | SMART 5    |         | IR         |         | Reference Oven |         |
|----------------------|------------|---------|------------|---------|------------|---------|----------------|---------|
|                      | Time (sec) | St. Dev | Time (sec) | St. Dev | Time (sec) | St. Dev | Time (hrs)     | St. Dev |
| Cream                | 135        | 0.05    | 210        | 0.05    | 600        | 0.1     | 3.5            | 0.03    |
| Sweet Condensed Milk | 150        | 0.11    | 270        | 0.18    | 720        | 0.25    | 4              | 0.06    |
| Cheese               | 120        | 0.12    | 240        | 0.15    | 1200       | 0.20    | 5              | 0.09    |
| Raw Meat             | 120        | 0.1     | 180        | 0.15    | 900        | 0.2     | 18             | 0.06    |
| Deli Meat            | 130        | 0.08    | 210        | 0.12    | 840        | 0.18    | 18             | 0.05    |
| Milk Powder          | 75         | 0.04    | 300        | 0.15    | 560        | 0.10    | 5              | 0.02    |
| Salad Dressing       | 110        | 0.02    | 200        | 0.03    | 900        | 0.03    | 5              | 0.01    |
| Latex                | 60         | 0.05    | 180        | 0.07    | 720        | 0.07    | 2.5            | 0.03    |
| Solvent Adhesive     | 180        | 0.12    | N/A        | N/A     | 480        | 0.22    | 2              | 0.1     |
| Coating              | 160        | 0.08    | 225        | 0.12    | 480        | 0.18    | 1              | 0.06    |

## THE SMART 6 CAN ANALYZE THESE AND MANY MORE...

Cod (Raw)  
 Salmon (Raw)  
 Heavy Cream  
 Milk (Past. Homog.)  
 Yogurt  
 Hot Dogs  
 Raw Pork  
 Chicken (Raw)  
 Deli Chicken (Cooked)

Turkey  
 Deli Ham  
 Ground Beef  
 Sausage  
 Bacon  
 Processed Cheese  
 Natural Cheese  
 Butter  
 Ice Cream Mix

Dressings  
 Mayonnaise  
 Condensed Skim Milk  
 Cream Cheese  
 Crackers  
 Chips  
 Noodles  
 Whey Powder

Dry Dog Food  
 Whole Milk Powder  
 Powdered Infant Formula  
 Cream Powder  
 Non-Fat Milk Powder

Wheat Flour  
 Poultry Feed  
 Processed Meat  
 Egg Whites  
 Duck  
 Lamb  
 Mixed Foods

Resin  
 Adhesives  
 Cosmetics  
 Detergents  
 Wastewater  
 Black Liquor  
 TiO<sub>2</sub>

# Ease of Use



## Individual User Logins

Password protected and CFR 21 Part 11 compliant



## 3 Simple Steps To Run Any Sample

1. Tare pads
2. Add Sample
3. Press Start



## Built-in Training Videos

Ensure your technicians run the SMART 6 correctly



## Method Import

Import from other users or CEM's Method Library from anywhere in the world



## Data Export

For easy data review and storage



# Safety

## Method Controls

- More controls over sample heating
- Ramps and stages to prevent any overheating or burning

## iPower Heating

- More efficient at heating samples than outdated technologies
- Chooses the best heating process for each specific sample

## Adjustable Fan Speeds

- Remove hazardous volatiles more quickly to create a safer heating environment
- Exhaust hose to safely move volatiles to an approved exhaust area

## Flame Detection

- Built-in IR sensor will detect sudden heat increases and immediately stop all heating if temperature gets too high
- Stops the fire before it even happens



# Upgrades/Options



## ORACLE™ Fat Analyzer

Adds AOAC approved methodology for fat analysis to the SMART 6. Get results in under 1 minute with no solvents or costly calibration maintenance.



## ProFat™ Meat Analyzer

Rapidly analyze raw meats and pre-blends for moisture, fat, and protein content in a single at-line system without the need for calibrations.



## Communicator Data Management

Automatically collect data, sort and analyze statistics based on dates, methods, and results, or chart results to easily view out-of-spec data.



## Bar Code Reader

Easy data tracking and sample identification with the optional barcode reader.



## Integrated Printer

Print results, method information, and system information with the optional thermal printer.

## In the box

- ✓ 200 glass fiber sample pads
- ✓ 2 Balance Pans
- ✓ 2 Balance Stems
- ✓ 1 Sample Press
- ✓ Monitor Kit
- ✓ 5 Pipettes
- ✓ 10 g Calibration Weight
- ✓ The Importance of Line Voltage booklet
- ✓ Setup Guide
- ✓ Power Cord (not pictured)

# CEM

3100 Smith Farm Road  
Matthews, NC 28104 USA

(704) 821-7015  
(800) 726-3331  
[www.cem.com](http://www.cem.com)

## EXPECT THE BEST

CEM will strive to meet or exceed our customers expectations through continuous improvement in pursuit of world class excellence. We will nurture the commitment to quality improvement in all aspects of our operations.

M.J. Collins  
President/CEO

# SMART<sub>6</sub>





Over 50,000  
systems sold  
worldwide



CEM has been an  
ISO-certified facility  
since 1994



All systems serviced  
& supported by  
experts with an  
average of 15 years  
of experience



CEM invests 11%  
of annual revenue  
into R&D, the result...  
11 R&D 100 awards



IQ/OQ/PQ Validation  
by certified CEM  
Technicians

“

*Our commitment to you  
doesn't end when your  
system is shipped;  
it begins.*

~Michael J. Collins  
President & CEO, CEM

”

**CEM Corporation:**  
PO Box 200  
Matthews, NC 28106  
United States

800-726-3331  
704-821-7015  
Fax: 704-821-7894  
info@cem.com  
www.cem.com

**France: CEM µWave S.A.S.**  
Immeuble Ariane  
Domaine Technologique de Saclay  
4, rue Rene' Razel, 91892 ORSAY Cedex  
33 (01) 69 35 57 80 • Fax: 33 (01) 60 19 64 91  
info.fr@cem.com

**Germany, Austria, Switzerland: CEM GmbH**  
Carl-Friedrich-Gauss-Str.9, 47475 Kamp-Lintfort  
(49) 2842-9644-0 • Fax: (49) 2842-9644-11  
info@cem.de

**Ireland: CEM Technology (Ireland) Ltd.**  
Sky Business Centre, 9a Plato Business Park,  
Damastown, Dublin 15  
+353 (0) 1 885 1752 • Fax: +353 (0) 1 885 1601  
info.ireland@cem.com

**Italy: CEM S.R.L.**  
Via Dell' Artigianato, 6/8  
24055 Cologno al Serio (Bg)  
(39) 35-896224 • Fax: (39) 35-891661  
info.srl@cem.com

**Japan: CEM Japan K.K.**  
2-18-10 Takanawa, Minato-ku, Tokyo 108-0074  
+81-3-5793-8542 • Fax: +81-3-5793-8543  
info@cemjapan.co.jp

**United Kingdom: CEM Microwave Technology Ltd.**  
2 Middle Slade, Buckingham Industrial Estate,  
Buckingham MK181WA  
(44) 1280-822873 • Fax: (44) 1280-822873  
info.uk@cem.com

6,227,041 6,320,170 6,521,876 6,268,570

Worldwide patents issued and pending

©2016 CEM Corporation

SMART 6™, Octawave™, iPower®, iDri™, Intelli-Temp™ are trademarks of CEM Corporation.

**CEM**