

# ELGA LabWater Product, Consumables & Accessory Guide



**This guide is not for use in USA**

## Deionisation Cylinders

Part Number	Description
<b>30 Litres per hour, 1-10MΩ-cm</b>	
B114/A	Elgastar B114 Deioniser Complete With Four C114 Cartridges
<b>Consumable</b>	
C114/4	Replacement Pack Of Four C114 Cartridges
<b>60 Litres per hour, 1-10MΩ-cm</b>	
VISION125	Wall-mounted deioniser
<b>Accessory and consumable</b>	
SP404	Dispense Gun and Line for Vision & TIU Unit
VC125	Replacement disposable cartridge
<b>60 Litres per hour, 1-10MΩ-cm</b>	
VISION250	Wall-mounted deioniser
<b>Accessory and consumable</b>	
SP404	Dispense Gun and Line for Vision & TIU Unit
VC250	Replacement disposable cartridge
<b>90 Litres per hour, 1-10MΩ-cm</b>	
M/M/DS/A	Elgastat Micromeg Deioniser complete with two MC:DS Cartridges
<b>Consumable</b>	
MC:DS/4	Replacement pack of four MC:DS Cartridges
<b>240 Litres per hour, 1-10MΩ-cm</b>	
B118/A	Deioniser complete with one C118 Cylinder
<b>240 Litres per hour, 1-50µS-cm</b>	
B119/A	Deioniser complete with one C119 Cylinder
<b>540 Litres per hour, 10-50µS-cm</b>	
B500/A	Deioniser complete with one C500 Cylinder
<b>540 Litres per hour, 1-10MΩ-cm</b>	
B510/A	Deioniser complete with one C510 Cylinder
<b>Consumable</b>	
MAXI/S	Maxipack S Refill, 62 litres of standard ion-exchange resin



This guide is not for use in USA

# PURELAB Prima

Part Number	Description
PR007BPM1	PURELAB Prima 7 with Raw Water Boost Pump
PR007XXM1	PURELAB Prima 7
PR015BPM1	PURELAB Prima 15 with Raw Water Boost Pump
PR015XXM1	PURELAB Prima 15
PR030BPM1	PURELAB Prima 30 with Raw Water Boost Pump
PR030XXM1	PURELAB Prima 30

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, Starter Pack of CT tablets, LA634 installation kit, manual

## Prima 7/15/30 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## Accessories for Prima 7/15/30

LA513	Installation Kit multi-purpose
LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LA582	Pre-filter assembly
LC123	5.0 Micron Briquetted Carbon Filter (10")
LA618	Printer Kit for Ultra and SBP range
LA610	Wall Mounting Kit SBP and Ultra
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA620	DV Modular Reservoir (35 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB Prima

Part Number	Description
PRI060XXXX11503-1	PURELAB Prima 60 (115V)
PRI060XXXX23003-1	PURELAB Prima 60 (230V)
PRI120XXXX11503-1	PURELAB Prima 120 (115V)
PRI120XXXX23003-1	PURELAB Prima 120 (230V)
PRI090XXXX11503-1	PURELAB Prima 90 (115V)
PRI090XXXX23003-1	PURELAB Prima 90 (230V)

Unit is supplied as standard with:

Fitted / Included: 1 x LC123, Appropriate part number and quantity RO units for unit size, LA545 installation kit, manual

## Prima 30/60/90/120 Consumables

LC123 5.0 Micron Briquetted Carbon Filter (10")

Recommended Change: Maximum 3 months

LC115 30L/Hr Reverse Osmosis Module

LC119 60L/Hr Reverse Osmosis Module

LC120 90L/Hr Reverse Osmosis Module

Recommended Change: all RO's typically 2-3 years

LC136M2 Composite Air Vent Filter

Recommended Change: 6 months

## Accessories for Prima 30/60/90/120

LA518 SP311 with 8/12mm John Guest Fittings

LA543 Soft Water Kit

LA593 Wall Mounting Kit (for unit)

LA545 Installation kit

FILT/JAC Backwashable Multi media filter (230V /50Hz)

FILT/JFE Backwashable Iron removal filter (230V /50Hz)

TEST30679 Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA586-1 25 litre Storage Tank

LA591 Wall-mounting Kit

LA587-1 40 litre Storage Tank

LA591 Wall-mounting Kit

LA590-1 75 litre Storage Tank

LA592 Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB Option-S

Part Number	Description
OS007BPM1	PURELAB Option S 7 with Raw Water Boost Pump
OS007XXM1	PURELAB Option S 7
OS015BPM1	PURELAB Option S 15 with Raw Water Boost Pump
OS015XXM1	PURELAB Option S 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC141, LA634 installation kit, manual

## Option-S 7/15 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC141	Twin-pack Cartridge
Recommended Change:	Maximum 6 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## Accessories for Option-S 7/15

LA513	Installation Kit multi-purpose
LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LC123	5.0 Micron Briquetted Carbon Filter (10")
LA582	Pre-filter assembly
LA604	Flow Kit-Option S 7 To 15 Conversion
LA610	Wall Mounting Kit SBP and Ultra
LA622	Wall Mounting Kit SBP (USA Version)
LA618	Printer Kit for Ultra and SBP range
LA609	Boost Pump Kit SBP
LA624	
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA620	DV Modular Reservoir (35 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB Option-S

Part Number	Description
OPT030XXXX11503-1	PURELAB Option S 30 (115V)
OPT030XXXX23003-1	PURELAB Option S 30 (230V)
OPT060XXXX11503-1	PURELAB Option S 60 (115V)
OPT060XXXX23003-1	PURELAB Option S 60 (230V)

Unit is supplied as standard with:

Fitted / Included: 1 x LC123, Appropriate part number and quantity RO units for unit size, 1 x LC114, 1 x LC105, bypass block, LA545 installation kit, manual

## Option-S 30/60 Consumables

LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	Maximum 3 months
LC115	30L/Hr Reverse Osmosis Module
LC119	60L/Hr Reverse Osmosis Module
Recommended Change:	Typically 2-3 years
LC114	Ion/Organic Removal Twin Pack
Recommended Change:	As required, maximum 6 months
LC105	UV Lamp
Recommended Change:	12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets

## Accessories for Option-S 30/60

LA518	SP311 with 8/12mm John Guest Fittings
LA543	Soft Water Kit
LA593	Wall Mounting Kit (for unit)
LA545	Installation kit
FILT/JAC	Backwashable Multi media filter (230V /50Hz)
FILT/JFE	Backwashable Iron removal filter (230V /50Hz)
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA586-1	25 litre Storage Tank
LA591	Wall-mounting Kit
LA587-1	40 litre Storage Tank
LA591	Wall-mounting Kit
LA590-1	75 litre Storage Tank
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB Pulse

Part Number	Description
PP010BPM1	PURELAB Pulse 1 with Raw Water Boost Pump
PP010XXM1	PURELAB Pulse 1 – (Up to 10L/hr)
PP020BPM1	PURELAB Pulse 2 with Raw Water Boost Pump
PP020XXM1	PURELAB Pulse 2 – (Up to 20L/hr)

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC178, 1 x LC105, 1 x LC218, Starter Pack of CT tablets & bypass block, LA636 installation kit, manual

## PURELAB Pulse Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC178	Conditioning Cartridge
Recommended Change:	6 months
LC105	UV Lamp
Recommended Change:	12 months
LC218	EDI Stack
Recommended Change:	4-5 Years
LC197	Point of Use Biofilter
Recommended Change:	6 months
LC145	Point of Use Filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## Accessories for PURELAB Pulse

LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LC123	5.0 Micron Briquetted Carbon Filter (10")
LA582	Pre-filter assembly
LA610	Wall Mounting Kit SBP and Ultra
LA513	Installation Kit multi-purpose
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA737	25ltr Docking vessel – DV25
LA591	Wall-mounting Kit
LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

## PURELAB Option-E

Part Number	Description
ELE025XXXX23003-1	PURELAB Option E 25 (230V /50Hz)
ELE050XXXX23003-1	PURELAB Option E 50 (230V /50Hz)
ELE080XXXX23003-1	PURELAB Option E 80 (230V /50Hz)

Unit is supplied as standard with:

Fitted / Included: 1 x LC123, Appropriate part number and quantity RO units for unit size, 1 x LC122, 1 x LC105, LA545 installation kit, manual

### Option-E 25/50/80 Consumables

LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	Maximum 3 months
LC115	30L/Hr Reverse Osmosis Module
LC119	60L/Hr Reverse Osmosis Module
Recommended Change:	Typically 2-3 years
LC122	Conditioning cartridge
Recommended Change:	3 months
LC105	UV Lamp
Recommended Change:	12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

### Accessories for Option-E 25/50/80

LA518	SP311 with 8/12mm John Guest Fittings
LA593	Wall Mounting Kit (for unit)
LA545	Installation kit
FILT/JAC	Backwashable Multi media filter (230V /50Hz)
FILT/JFE	Backwashable Iron removal filter (230V /50Hz)
TEST30679	Chlorine Test Kit

### Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA586-1	25 litre Storage Tank
LA591	Wall-mounting Kit
LA587-1	40 litre Storage Tank
LA591	Wall-mounting Kit
LA590-1	75 litre Storage Tank
LA592	Wall Mounting Kit



**This guide is not for use in USA**



# PURELAB Option-R

Part Number	Description
OR007BPM1	PURELAB Option R 7 with Raw Water Boost Pump
OR007XXM1	PURELAB Option R 7
OR015BPM1	PURELAB Option R 15 with Raw Water Boost Pump
OR015XXM1	PURELAB Option R 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC105, 1 x LC141, Starter Pack of CT tablets & sanitisation block, LA635 installation kit, manual

## Option-R 7/15 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC105	UV Lamp
Recommended Change:	12 months
LC141	Twin-pack Cartridge
Recommended Change:	Maximum 6 months
LC197	Point of Use Biofilter
Recommended Change:	6 months
LC145	Point of Use Filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## Option-R 7/15 Accessories

LA513	Installation Kit multi-purpose
LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LC123	5.0 Micron Briquetted Carbon Filter (10")
LA582	Pre-filter assembly
LA605	Flow Kit-Option R 7 To 15 Conversion
LA610	Wall Mounting Kit SBP and Ultra
LA618	Printer Kit for Ultra and SBP range
LA625	Remote Control Station
LA629	Remote Dispense Gun
LA643	Remote Dispense Station (Including Dispense Gun)
LA521	Remote Dispense Valve
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

## PURELAB Option-R

Part Number	Description
OPT030PLXX11503-1	PURELAB Option R 30 (115V)
OPT030PLXX23003-1	PURELAB Option R 30 (230V)
OPT060PLXX11503-1	PURELAB Option R 60 (115V)
OPT060PLXX23003-1	PURELAB Option R 60 (230V)

Unit is supplied as standard with:

Fitted / Included: 1 x LC123, Appropriate part number and quantity RO units for unit size, 1 x LC114, 1 x LC105, bypass block, LA545 installation kit, manual

### Option-R 30/60 Consumables

LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	Maximum 3 months
LC115	30L/Hr Reverse Osmosis Module
LC119	60L/Hr Reverse Osmosis Module
Recommended Change:	Typically 2-3 years
LC114	Ion/Organic Removal Twin Pack
Recommended Change:	As required, maximum 6 months
LC105	UV Lamp
Recommended Change:	12 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets

### Accessories for Option-R 30/60

LA545	Installation kit
LA518	SP311 with 8/12mm John Guest Fittings
LA543	Soft Water Kit
LA593	Wall Mounting Kit (for unit)
FILT/JAC	Backwashable Multi media filter (230V /50Hz)
FILT/JFE	Backwashable Iron removal filter (230V /50Hz)
TEST30679	Chlorine Test Kit

### Tanks

LA586-1	25 litre Storage Tank
LA591	Wall-mounting Kit
LA587-1	40 litre Storage Tank
LA591	Wall-mounting Kit
LA590-1	75 litre Storage Tank
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB 3000

Part Number	Description
PB3060XXM1-115	PURELAB 3000 60LPH 115V
PB3060XXM1-230	PURELAB 3000 60LPH 230V
PB3120XXM1-115	PURELAB 3000 120LPH 115V
PB3120XXM1-230	PURELAB 3000 120LPH 230V

Unit is supplied as standard with:

Fitted / Included: 1 x LC202, Appropriate RO (LC180 or LC193), Quick Start reference and operator manuals, user recognition key set

## PURELAB 3000 Starter Kits

LA721	Starter Kit 1
Contains	1 x LC175, LA720 installation kit
LA722	Starter Kit 2
Contains	1 x LC177, LA720 installation kit

## PURELAB 3000 Consumables

LC175	Protek L1 Pretreatment Cartridge
Recommended Change:	6 months
LC177	Protek L2 Pretreatment Cartridge
Recommended Change:	12 months
LC180	60L/Hr Reverse Osmosis Cartridge
LC193	120 l/hr Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 Years
LC181	E - Cartridge
Recommended Change:	36 months
LC187	Purification Pack
Recommended Change:	3-12 months
LC202	Vent Filter
Recommended Change:	6 months
CT3	Pack of 12 Effervescent Sanitisation Tablets
CT3-50	Pack of 50 Effervescent Sanitisation Tablets
CT3-300	6 Packs of 50 Effervescent Sanitisation Tablets
Recommended Change:	Typical Usage is 2 Tablets per month

## PURELAB 3000 Accessories

LA701	E-Cartridge Upgrade Kit complete with LC181
LA703	60LPH to 120LPH Upgrade Kit
LA704	Castors
TEST30679	Chlorine Test Kit



**This guide is not for use in USA**

# PURELAB 7000

Part Number	Description
PB7060XXM1-115	PURELAB 7000 60LPH 115V
PB7060XXM1-230	PURELAB 7000 60LPH 230V
PB7120XXM1-115	PURELAB 7000 120LPH 115V
PB7120XXM1-230	PURELAB 7000 120LPH 230V

Unit is supplied as standard with:

Fitted / Included: 1 x LC136M2, Appropriate RO (LC180 or LC193), 1 x LC118, 1 x LC109, Quick Start reference and operator manuals, user recognition key set

## PURELAB 7000 Starter Kits

LA729	Starter Kit 1
Contains	1 x LC175, 1 x LC187, LA720 installation kit
LA730	Starter Kit 2
Contains	1 x LC177, 1 x LC187, LA720 installation kit

## PURELAB 7000 Consumables

LC175	Protek L1 Pretreatment Cartridge
Recommended Change:	6 months
LC177	Protek L2 Pretreatment Cartridge
Recommended Change:	12 months
LC180	60L/Hr Reverse Osmosis Cartridge
LC193	120 l/hr Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 Years
LC181	E - Cartridge
Recommended Change:	36 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months
LC118	UV Lamp
Recommended Change:	12 months
LC187	Purification Pack
Recommended Change:	3-12 months
LC109	0.05 Micron ultra-microfilter (UMF)
Recommended Change:	6 months
CT3	Pack of 12 Effervescent Sanitisation Tablets
CT3-50	Pack of 50 Effervescent Sanitisation Tablets
CT3-300	6 Packs of 50 Effervescent Sanitisation Tablets
Recommended Change:	Typical Usage is 2 Tablets per month

## PURELAB 7000 Accessories

LA723	High Flow Upgrade Kit complete with LC187
LA701	E-Cartridge Upgrade Kit complete with LC181
LA703	60LPH to 120LPH Upgrade Kit
LA704	Castors
LA631	40 Litre Tank complete with Level Control Assembly 230V 50Hz
LA633	40 Litre Tank complete with Level Control Assembly 115V 60Hz
LA643	Remote Dispense Station (Including Dispense Gun)
TEST30679	Chlorine Test Kit



**This guide is not for use in USA**

## PURELAB Ultra

Part Number	Description
ULRSXANM2	PURELAB Ultra Analytic under bench system with remote control station and dispense gun.
ULRSXGEM2	PURELAB Ultra Genetic under bench system with remote control station and dispense gun.
ULXXXANM2	PURELAB Ultra Analytic ultra-pure polishing system.
ULXXXBSM2	PURELAB Ultra Bioscience ultra-pure polishing system.
ULXXXGEM2	PURELAB Ultra Genetic ultra-pure polishing system.
ULXXXIOM2	PURELAB Ultra Ionic ultra-pure polishing system.
ULXXXSCM2	PURELAB Ultra Scientific ultra-pure polishing system.

Unit is supplied as standard with:

Fitted / Included: Appropriate consumables (EXCEPT Purification Packs), Bypass block and starter pack of sanitisation tablets, LA642 installation kit including inlet strainer, manual

### Range of Purification Packs

LC182	Labpure S1 Purification Cartridge - RO Feed
LC183	Labpure S2 Purification Cartridge - SDI Feed
LC184	Labpure S3 Purification Cartridge - Low Ionic
LC185	Labpure S4 Purification Cartridge - Low TOC
<i>Recommended Change: Maximum 6 months</i>	

### Consumables - in Addition to Purification Cartridges

#### PURELAB Ultra Scientific Consumables

LC134	Point of Use filter
<i>Recommended Change: 6 months</i>	
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change: Use 1 tablet, once a month. Pack size therefore 12 months</i>	

#### PURELAB Ultra Ionic Consumables

LC105	UV Lamp
<i>Recommended Change: 12 months</i>	
LC109	0.05 Micron ultra-microfilter (UMF)
<i>Recommended Change: 6 months</i>	
LC197	Point of Use Biofilter
<i>Recommended Change: 6 months</i>	
LC134	Point of Use filter
<i>Recommended Change: 6 months</i>	
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change: Use 1 tablet, once a month. Pack size therefore 12 months</i>	

#### PURELAB Ultra Analytic Consumables

LC118	UV Lamp
<i>Recommended Change: 12 months</i>	
LC109	0.05 Micron ultra-microfilter (UMF)
<i>Recommended Change: 6 months</i>	
LC134	Point of Use filter
<i>Recommended Change: 6 months</i>	
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change: Use 1 tablet, once a month. Pack size therefore 12 months</i>	



**This guide is not for use in USA**

## PURELAB Ultra

### Part Number Description

#### PURELAB Ultra Bioscience Consumables

LC151	Endoguard 5000 UF Membrane
<i>Recommended Change:</i>	<i>24 months</i>
LC134	Point of Use filter
<i>Recommended Change:</i>	<i>6 months</i>
LC197	Point of Use Biofilter
<i>Recommended Change:</i>	<i>6 months</i>
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change:</i>	<i>Use 1 tablet, once a month. Pack size therefore 12 months</i>

#### PURELAB Ultra Genetic Consumables

LC118	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC151	Endoguard 5000 UF Membrane
<i>Recommended Change:</i>	<i>24 months</i>
LC134	Point of Use filter
<i>Recommended Change:</i>	<i>6 months</i>
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change:</i>	<i>Use 1 tablet, once a month. Pack size therefore 12 months</i>

#### Tanks *(supplied as standard with 1 x LC136M2 composite vent filter)*

LA620	DV Modular Reservoir (35 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit
LC136M2	Composite Air Vent Filter
<i>Recommended Change:</i>	<i>6 months</i>

#### PURELAB Ultra Accessories

LA610	Wall Mounting Kit SBP and Ultra
LA622	Wall Mounting Kit SBP (USA Version)
LA646	Open Circuit Ringmain Kit
LA647	Closed Circuit Ringmain Kit
LA618	Printer Kit for Ultra and SBP range
LA575	Pressure regulator with relief line
LA652	Inlet Pressure Regulator
LA642	Installation Kit
LA644	Integral Dispense Gun
LA643	Remote Dispense Station (Including Dispense Gun)
LA645	Remote Control Station
LA698	Sanitisation by pass block
TEST30679	Chlorine Test Kit



**This guide is not for use in USA**

## PURELAB Classic

Part Number	Description
CLXXUVM2	PURELAB Classic UV UF
CLXXDIM2	PURELAB Classic DI
CLXXUVM2	PURELAB Classic UF
CLXXUVM2	PURELAB Classic UV

Unit is supplied as standard with:

Fitted / Included: Appropriate consumables (EXCEPT Purification Packs), Bypass block and starter pack of sanitisation tablets, LA642 installation kit including inlet strainer, manual

### PURELAB Classic DI Consumables

LC186	Labpure S5 Purification Cartridge
Recommended Change:	Maximum 12 months
LC134	Point of Use filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months

### PURELAB Classic UF Consumables

LC186	Labpure S5 Purification Cartridge
Recommended Change:	Maximum 12 months
LC169	UF Membrane
Recommended Change:	24 months
LC134	Point of Use filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months

### PURELAB Classic UV Consumables

LC170	UV Lamp
Recommended Change:	12 months
LC186	Labpure S5 Purification Cartridge
Recommended Change:	Maximum 12 months
LC134	Point of Use filter
Recommended Change:	6 months
LC197	Point of Use Biofilter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months

### PURELAB Classic UVF Consumables

LC170	UV Lamp
Recommended Change:	12 months
LC186	Labpure S5 Purification Cartridge
Recommended Change:	Maximum 12 months
LC169	UF Membrane
Recommended Change:	24 months
LC134	Point of Use filter
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, once a month. Pack size therefore 12 months



**This guide is not for use in USA**

## PURELAB Classic

Part Number	Description
<b>Tanks</b>	
<i>(supplied as standard with 1 x LC136M2 composite vent filter)</i>	
LA620	DV Modular Reservoir (35 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit
LC136M2	Composite Air Vent Filter

Recommended Change: 6 months

### PURELAB Classic Accessories

LA610	Wall Mounting Kit SBP and Ultra
LA618	Printer Kit for Ultra and SBP range
LA575	Pressure regulator with relief line
LA652	Inlet Pressure Regulator
LA642	Installation Kit
LA646	Open Circuit Ringmain Kit
LA647	Closed Circuit Ringmain Kit
LA644	Integral Dispense Gun
LA645	Remote Control Station
LA698	Sanitisation by pass block



**This guide is not for use in USA**



# PURELAB Option-Q

Part Number	Description
OQ007BPM1	PURELAB Option Q 7 with Raw Water Boost Pump
OQ007XXM1	PURELAB Option Q 7
OQ015BPM1	PURELAB Option Q 15 with Raw Water Boost Pump
OQ015XXM1	PURELAB Option Q 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC118, 1 x LC163, Starter Pack of CT tablets & sanitisation block, LA635 installation kit, manual

## PURELAB Option-Q 7/15 Consumables

LC140	Pre-treatment Cartridge
<i>Recommended Change:</i>	<i>Maximum 6 months</i>
LC143	Reverse Osmosis Cartridge
<i>Recommended Change:</i>	<i>Typically 2-3 years</i>
LC118	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC163	Twin Pack Purification Cartridge
<i>Recommended Change:</i>	<i>6 months</i>
LC145	Point of Use Filter
<i>Recommended Change:</i>	<i>6 months</i>
LC197	Point of Use Biofilter
<i>Recommended Change:</i>	<i>6 months</i>
CT1	Pack Of 12 Disinfection Tablets
<i>Recommended Change:</i>	<i>Use 1 tablet, once a month. Pack size therefore 12 months</i>
LC136M2	Composite Air Vent Filter
<i>Recommended Change:</i>	<i>6 months</i>

## PURELAB Option-Q Accessories

LC123	5.0 Micron Briquetted Carbon Filter (10")
LA513	Installation Kit multi-purpose
LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LA582	Pre-filter assembly
LA605	Flow Kit-Option R 7 To 15 Conversion
LA610	Wall Mounting Kit SBP and Ultra
LA622	Wall Mounting Kit SBP (USA Version)
LA618	Printer Kit for Ultra and SBP range
LA675	Remote Control Station
LA643	Remote Dispense Station (Including Dispense Gun)
LA521	Remote Dispense Valve
TEST30679	Chlorine Test Kit

## Tanks *(supplied as standard with 1 x LC136M2 composite vent filter)*

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# PURELAB flex

Part Number	Description
PF1XXXXM1	PURELAB flex (optional Purification Pack)
PF2XXXXM1	PURELAB flex 2 (Purification Pack + UV + TOC)
PF3XXXXM1	PURELAB flex 3
PF4XXXXM1	PURELAB flex 4
PF5XXXXM1	PURELAB flex 5
PF6XXXXM1	PURELAB flex 6

Unit is supplied as standard with:

All units: Quick Reference Guide & user manual, certificate of conformity, bypass pack.

Flex 2: 1 x LC210. Flex 3 & 5: 1 x LC217, 1 x LC210, 1 x LC216. Flex 4 & 6: 1 x LC210, 1 x LC216.

## Consumables for flex 1 and 2

LC208	Purification Pack
Recommended Change:	12 months
LC209	Sanitisation Pack
Recommended Change:	Single use sanitisation pack
LC145	Point of Use Filter
Recommended Change:	3 months
LC197	Point of Use Biofilter
Recommended Change:	3 months

## For flex PF2XXXXM1 only

LC210	UV Lamp
Recommended Change:	12-18 months

## Accessories:

LA731	Pre-Conditioning Kit
LC211	Pre-Conditioning Cartridge
LA732	Foot Switch
LA734	Leak Detector/Controller
LA652	Inlet Pressure Regulator
LA728	Bypass Pack - Flex

## Consumables for flex 3 - 6

LC214	Purification Pack
Recommended Change:	6 months
LC209	Sanitisation Pack
Recommended Change:	Single use sanitisation pack
LC210	UV Lamp
Recommended Change:	12-18 months
LC216	Composite Vent Filter
Recommended Change:	12 months
LC145	Point of Use Filter
Recommended Change:	3 months
LC197	Point of Use Biofilter
Recommended Change:	3 months

## For flex 3 & 5 only

LC217	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years

## Accessories

LA732	Foot Switch
LA734	Leak Detector/Controller
LA735	Wall Mounting Kit
LA652	Inlet Pressure Regulator
LA728	Bypass Pack - Flex



**This guide is not for use in USA**

## PURELAB Chorus

### Part Number Description

#### PURELAB Chorus 3

RO310XXM1	PURELAB Chorus 3 RO 10l/Hr
RO310BPM1	PURELAB Chorus 3 RO 10l/Hr With Boost Pump
RO320XXM1	PURELAB Chorus 3 RO 20l/Hr
RO320BPM1	PURELAB Chorus 3 RO 20l/Hr With Boost Pump
RO330XXM1	PURELAB Chorus 3 RO 30l/Hr
RO330BPM1	PURELAB Chorus 3 RO 30l/Hr With Boost Pump

Unit is supplied as standard with:

Fitted / Included: Appropriate quantity LC240 RO, 1 x LC241, GUID39864 Light Guide, GUID40004 Quick Reference Guide, INST40008 Quick Start Guide, LA762 basic installation kit

#### PURELAB Chorus 3 Consumables

LC241	Pre-Treatment Pack
Recommended Change:	6 months
LC240	10L/Hr Reverse Osmosis Module
Recommended Change:	Typically 2-3 Years

#### PURELAB Chorus 2

PC2RODIM1	PURELAB Chorus 2 RO/DI System
-----------	-------------------------------

Unit is supplied as standard with:

Fitted / Included: Appropriate quantity LC240 RO, 1 x LC241, 1 x LC234, MANU40003, Operator Manual, GUID39864 Light Guide, GUID40005 Quick Reference Guide, INST40009 Quick Start Guide, LA762 basic installation kit, LC233 bypass block

#### PURELAB Chorus 2 Consumables

LC241	Pre-Treatment Pack
Recommended Change:	6 months
LC240	10L/Hr Reverse Osmosis Module
Recommended Change:	Typically 2-3 Years
LC234	Purification Pack - RO/DI
Recommended Change:	12 months

#### PURELAB Chorus Reservoirs

LA757	15 Litre Reservoir
LA758	30 Litre Reservoir
LA759	60 Litre Reservoir

Each reservoir supplied with LA773 Installation Kit & LC216 Composite Vent Filter

#### PURELAB Chorus Reservoir Consumable

LC216	Composite Vent Filter
Recommended Change:	6 months

#### PURELAB Chorus Reservoir Accessories

LA770	Wall Mounting Kit 15/30 L Reservoir
LA771	Wall Mounting Kit 60/100 L Reservoir

#### PURELAB Chorus 1

PC1ANRXM1	PURELAB Chorus 1 Analytic Research Ultrapure System
PC1LSCXM1	PURELAB Chorus 1 Life Science Ultrapure System
PC1GSCXM1	PURELAB Chorus 1 General Science Ultrapure System

Unit is supplied as standard with:

Fitted / Included: Appropriate consumables (including 2 x LC232 - see list below), MANU39998/MANU39997/MANU40001 Operator Manual, GUID39864 Light Guide, GUID40005 Quick Reference Guide, INST40012 Quick Start Guide, LA762 basic installation kit, LC233 x 2 bypass block



**This guide is not for use in USA**

## PURELAB Chorus

### Part Number Description

#### Purification Pack options

LC232	Purification Pack - RO Feed
LC244	Purification Pack - SDI Feed
LC245	Purification Pack - Low Ionic
LC246	Purification Pack - Low TOC
<i>Recommended Change:</i>	<i>12 months</i>
LC236	Sanitisation Pack (Closed Loop)
<i>Recommended Change:</i>	<i>Single Use</i>

**PURELAB Chorus General Science** only uses purification pack

#### PURELAB Chorus Analytic Consumables (additional to purification packs)

LC210	UV Lamp
<i>Recommended Change:</i>	<i>12-18 months</i>
LC109	0.05 Micron ultra-microfilter (UMF)
<i>Recommended Change:</i>	<i>12 months</i>

#### PURELAB Chorus Life Science Consumables (additional to purification packs)

LC210	UV Lamp
<i>Recommended Change:</i>	<i>12-18 months</i>
LC151	Endoguard 5000 UF Membrane
<i>Recommended Change:</i>	<i>12 months</i>

#### PURELAB Chorus Halo Dispensers

LA754	Halo Dispenser
LA755	Advanced Halo Dispenser
LA756	Flexible Dispenser

*Each supplied with GUID40007 or GUID40006 Quick Reference Guide, INST40011 or INST40013 Quick Start Guide, LA774 Installation Kit (plus for LA756: MANU40002 Operator manual)*

#### PURELAB Chorus Halo Accessories

LA768	Remote Mounting Kit For Halo Dispenser (All Models)
-------	---

#### PURELAB Chorus Halo Dispenser Consumables

LC145	Point of Use Filter
<i>Recommended Change:</i>	<i>3 months</i>
LC197	Point of Use Biofilter
<i>Recommended Change:</i>	<i>3 months</i>

#### PURELAB Chorus Installation Kits & Accessories

CEPWGU303754	PURELAB Chorus Power Supply Desktop 160W
LA762	Basic Installation Kit (Tubing & Bends)
LA773	Reservoir (All Sizes) Installation Kit
LA774	Halo Installation Kit
LA777	Integrated Installation Kit
LA778	Side By Side Installation Kit
LA779	Remote Mount Installation Kit
LA776	Tool Kit
LA512	Pressure regulator
LA652	Inlet Pressure Regulator
LA575	Pressure regulator with relief line
LC233	Bypass Block
LA775	Degassing kit for Chorus 3 (RO) and 2 (RO/DI)
LC181	E - Cartridge
<i>Recommended Change:</i>	<i>3 Years</i>
LA732	Foot Switch



**This guide is not for use in USA**

# MEDICA-S

Part Number	Description
ME007BPM1	MEDICA 7 with Raw Water Boost Pump
ME007PK01	MEDICA 7 BP Pack complete with LA621
ME007XXM1	MEDICA 7
ME015BPM1	MEDICA 15 with Raw Water Boost Pump
ME015XXM1	MEDICA 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC105, 1 x LC141, 1 x LC125, Starter Pack of CT tablets & sanitisation block, LA637 installation kit, manual. ME007PK01 also supplied with LA621 & LC136M2

## MEDICA 7/15 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC105	UV Lamp
Recommended Change:	12 months
LC141	Twin-pack Cartridge
Recommended Change:	Maximum 6 months
LC125	0.2 Micron in line capsule filter with adaptors
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, every 1-6 months, dependant on local conditions.
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## MEDICA 7/15 Accessories:

LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	6 months
LA582	Pre-filter assembly
LA610	Wall Mounting Kit SBP and Ultra
LA618	Printer Kit for Ultra and SBP range
LA606	Flow Kit-MEDICA 7 To 15 Conversion
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

# MEDICA-R

Part Number	Description
MR007BPM1	MEDICA R 7 with Raw Water Boost Pump
MR007PK02	MEDICA R 7 BP Pack complete with LA621
MR007XXM1	MEDICA R 7
MR015BPM1	MEDICA R 15 with Raw Water Boost Pump
MR015PK01	MEDICA R 15 Pack complete with LA621
MR015PK02	MEDICA R 15 BP Pack complete with LA621
MR015XXM1	MEDICA R 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC105, 1 x LC141, 1 x LC125, Starter Pack of CT tablets & sanitisation block, LA637 installation kit, manual. MR007PK01, MR015PK01, MR015PK02 also supplied with LA621, LA512, LA652 & LC136M2

## MEDICA-R 7/15 Consumables:

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC105	UV Lamp
Recommended Change:	12 months
LC141	Twin-pack Cartridge
Recommended Change:	Maximum 6 months
LC125	0.2 Micron in line capsule filter with adaptors
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, every 1-6 months, dependant on local conditions.
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## MEDICA-R 7/15 Accessories:

LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LA582	Pre-filter assembly
LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	6 months
LA610	Wall Mounting Kit SBP and Ultra
LA618	Printer Kit for Ultra and SBP range
LA606	Flow Kit-MEDICA 7 To 15 Conversion
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

## Tanks *(supplied as standard with 1 x LC136M2 composite vent filter)*

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**

## MEDICA-D

Part Number	Description
MD007BPM2	MEDICA D 7 complete with Raw Water Boost Pump
MD007XXM2	MEDICA D 7
MD015BPM2	MEDICA D 15 complete with Raw Water Boost Pump
MD015XXM2	MEDICA D 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC105, 1 x LC141, 1 x LC125, 1 x LC154, Starter Pack of CT tablets & sanitisation block, LA637 installation kit, manual

### MEDICA-D 7/15 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	Maximum 6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC141	Twin-pack Cartridge
Recommended Change:	Maximum 6 months
LC154	Degassing Membrane
Recommended Change:	Typically 2-3 years
LC105	UV Lamp
Recommended Change:	12 months
LC125	0.2 Micron in line capsule filter with adaptors
Recommended Change:	6 months
CT1	Pack Of 12 Disinfection Tablets
Recommended Change:	Use 1 tablet, every 1-6 months, dependant on local conditions.
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

### MEDICA-D 7/15 Accessories:

LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LA582	Pre-filter assembly
LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	6 months
LA610	Wall Mounting Kit SBP and Ultra
LA622	Wall Mounting Kit SBP (USA Version)
LA618	Printer Kit for Ultra and SBP range
LA606	Flow Kit-MEDICA 7 To 15 Conversion
LA609	Boost Pump Kit SBP
TEST30679	Chlorine Test Kit

### Tanks (supplied as standard with 1 x LC136M2 composite vent filter)

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit
LA611	25 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA612	40 Litre Tank with 5 level Switches
LA591	Wall-mounting Kit
LA613	75 Litre Tank with 5 level Switches
LA592	Wall Mounting Kit



**This guide is not for use in USA**



# MEDICA PRO

Part Number	Description
MP030RBM1-115	MEDICA Pro R 30 (115V)
MP030RBM1-230	MEDICA Pro R 30 (230V)
MP060RBM1-115	MEDICA Pro R 60 (115V)
MP060RBM1-230	MEDICA Pro R 60 (230V)
MP120RBM1-115	MEDICA Pro R 120 (115V)
MP120RBM1-230	MEDICA Pro R 120 (230V)

*Unit is supplied as standard with:*

*Fitted / Included: 1 x LC105, 1 x LC136M2, 1 x LC109, Appropriate RO type and Quantity for unit, Quick Start reference and operator manuals, user recognition key set*

MP030REM1-115	MEDICA Pro RE 30 (115V)
MP030REM1-230	MEDICA Pro RE 30 (230V)
MP060REM1-115	MEDICA Pro RE 60 (115V)
MP060REM1-230	MEDICA Pro RE 60 (230V)
MP120REM1-115	MEDICA Pro RE 120 (115V)
MP120REM1-230	MEDICA Pro RE 120 (230V)

*Unit is supplied as standard with:*

*Fitted / Included: 1 x LC105, 1 x LC136M2, 1 x LC109, 1 x LC181, Appropriate RO type and Quantity for unit, Quick Start reference and operator manuals, user recognition key set*

MPXXXLPM1	MEDICA Pro LPS
-----------	----------------

*Unit is supplied as standard with:*

*Fitted / Included: 1 x LC105, 1 x LC109, 1 x LC136M2, Appropriate RO type and Quantity for unit, Quick Start reference and operator manuals, user recognition key set*

MP060RXM1-115	MEDICA Pro EDI 60 (115V)
MP060RXM1-230	MEDICA Pro EDI 60 (230V)
MP120RXM1-115	MEDICA Pro EDI 120 (115V)
MP120RXM1-230	MEDICA Pro EDI 120 (230V)

*Unit is supplied as standard with:*

*Fitted / Included: 1 x LC105, 1 x LC136M2, Appropriate RO type and Quantity for unit, 1 x LC220 ADEPT Stack, 1 x LC181, Quick Start reference and operator manuals, user recognition key set*

## MEDICA Pro Starter Kits

LC188	Starter Kit 1: LA699, 1 LC175,1 LC174
LC189	Starter Kit 2: LA699, 1 LC177,1 LC174
LC190	Starter Kit 3 (for LPS only): LA699,1 LC192,1 LC174
LC198	Starter Kit 4: LA699, 2 LC174, 1 LC175, 1 CT3
LC199	Starter Kit 5: LA699, 2 LC175, 1 LC177, 1 CT3



**This guide is not for use in USA**



# MEDICA PRO

## Consumables

LC179	30L/Hr Reverse Osmosis Cartridge
LC180	60L/Hr Reverse Osmosis Cartridge
LC193	120 l/hr Reverse Osmosis Cartridge
<i>Recommended Change:</i>	<i>Typically 3 Years</i>
LC105	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC109	0.05 Micron ultra-microfilter (UMF)
<i>Recommended Change:</i>	<i>6 months</i>
LC174	Medpure L1 DI Pack
<i>Recommended Change:</i>	<i>3-12 months</i>
LC175	Protek L1 Pretreatment Cartridge
<i>Recommended Change:</i>	<i>6 months</i>
LC177	Protek L2 Pretreatment Cartridge
<i>Recommended Change:</i>	<i>12 months</i>
LC181	E - Cartridge
<i>Recommended Change:</i>	<i>3 Years</i>
CT3	Pack of 12 Effervescent Sanitisation Tablets
CT3-50	Pack of 50 Effervescent Sanitisation Tablets
CT3-300	6 Packs of 50 Effervescent Sanitisation Tablets
<i>Typical usage:</i>	<i>Use up to 6 tablets, every 1-6 months, dependant on local conditions.</i>
LC136M2	Composite Air Vent Filter
<i>Recommended Change:</i>	<i>12 months</i>

## Consumables for EDI units only

LC206	Conditioning Cartridge
<i>Recommended Change:</i>	<i>2-6 months</i>
LC207	Deionisation Cartridge
<i>Recommended Change:</i>	<i>3-12 months</i>
LC220	ADEPT 5 Stack
<i>Recommended Change:</i>	<i>Up to 5 years</i>

## Accessories

LA700	High Flow Upgrade Kit complete with LC174
LA701	E-Cartridge Upgrade Kit complete with LC181
LA702	30LPH to 60LPH Upgrade Kit
LA703	60LPH to 120LPH Upgrade Kit
LA704	Castors
LA643	Remote Dispense Station (Including Dispense Gun)



**This guide is not for use in USA**

## MEDICA R200

Part Number	Description
MR200DSM1-230	MEDICA R 200 (230V 50Hz) Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration
MR200DSM1-115	MEDICA R 200 (115V 60Hz) Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration

Unit is supplied as standard with:

Fitted / Included: 1 x LC156, 1 x LC159, 1 x LC158, 1 x LC160, 1 x LC161, Operator manual, PASSkey set, 2 x Filter spanners and 2 x sanitisation tubes

### Consumables

LC156	Composite Vent Filter
<i>Recommended Change:</i>	<i>6 months</i>
LC158	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC159	Pre-treatment (pump protection) Filter
<i>Recommended Change:</i>	<i>6 months</i>
LC160	0.2 micron filter
<i>Recommended Change:</i>	<i>6 months</i>
LC161	Reverse Osmosis Module
<i>Recommended Change:</i>	<i>Typically 3 Years</i>

### Accessories

LA665	Remote Display
LA667	Replacement Pass Key (Blue Operator Key)
LA668	Replacement Pass Key (Green Sanitisation Key)
LA679	High recovery (75%) Reverse Osmosis kit

### Aquadem Cylinder (example)

PRDISI0067	E 1000 Cylinder, with 94 litres of MB20 resins, fitted with High Flow Head
VMGUVA200684	High flow top cap with 3/4" Male Connectors
FPPIPV203045	Tubing set for use with High Flow Top Cap. 3/4" female Connectors



**This guide is not for use in USA**

# BIOPURE

Part Number	Description
MB007BPM1	BIOPURE 7 with Raw Water Boost Pump
MB007PK01	BIOPURE 7 BP Pack complete with LA621
MB007XXM1	BIOPURE 7
MB015BPM1	BIOPURE 15 with Raw Water Boost Pump
MB015PK01	BIOPURE 15 BP Pack complete with LA621
MB015XXM1	BIOPURE 15

Unit is supplied as standard with:

Fitted / Included: 1 x LC140, Appropriate quantity LC143 for unit size, 1 x LC105, 1 x LC191, 1 x LC151, 1 x CT3 tablets & sanitisation block & silicon tubing, LA717 installation kit, manual

## MEDICA Biopure 7/15 Consumables

LC140	Pre-treatment Cartridge
Recommended Change:	6 months
LC143	Reverse Osmosis Cartridge
Recommended Change:	Typically 2-3 years
LC105	UV Lamp
Recommended Change:	12 months
LC191	Optipure Purification Cartridge S1
Recommended Change:	Maximum 6 months
LC151	Endoguard 5000 UF Membrane
Recommended Change:	12 months
CT3	Pack of 12 Effervescent Sanitisation Tablets
Recommended Change:	Use 2 tablet each month. Pack size therefore 6 months
LC136M2	Composite Air Vent Filter
Recommended Change:	6 months

## MEDICA Biopure 7/15 Accessories

LA513	Installation Kit multi-purpose
LA506	Installation Kit with saddle valve for 15mm pipe
LA512	Pressure regulator
LA518	SP311 with 8/12mm John Guest Fittings
LC123	5.0 Micron Briquetted Carbon Filter (10")
Recommended Change:	6 months
LA591	Wall-mounting Kit
LA605	Flow Kit-Option R 7 To 15 Conversion
LA609	Boost Pump Kit SBP
LA521	Remote Dispense Valve
LA618	Printer Kit for Ultra and SBP range
TEST30679	Chlorine Test Kit
LA643	Remote Dispense Station (Including Dispense Gun)
LA714	Hygienic Sample Point (8mm & 10mm)
LA582	Pre-filter assembly

## Tanks

(supplied as standard with 1 x LC136M2 composite vent filter)

LA621	DV Modular Reservoir (25 litres)
LA591	Wall-mounting Kit



**This guide is not for use in USA**

# BIOPURE

Part Number	Description
MB060XXM1-115	BIOPURE 60 (115V)
MB060XXM1-230	BIOPURE 60 (230V)
MB120XXM1-115	BIOPURE 120 (115V)
MB120XXM1-230	BIOPURE 120 (230V)

Unit is supplied as standard with:

Fitted / Included: 1 x LC105, 1 x LC202, Appropriate RO type for unit, sanitisation bypass block, Quick Start reference guide, reference card and operator manual, user recognition key set

## MEDICA Biopure 60/120 Consumables

LC175	Protex L1 Pretreatment Cartridge
<i>Recommended Change:</i>	<i>6 months</i>
LC177	Protex L2 Pretreatment Cartridge
<i>Recommended Change:</i>	<i>12 months</i>
LC180	60L/Hr Reverse Osmosis Cartridge
LC193	120 l/hr Reverse Osmosis Cartridge
<i>Recommended Change:</i>	<i>Typically 3 years</i>
LC105	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC181	E - Cartridge
<i>Recommended Change:</i>	<i>3 years</i>
LC194	UF Filter
<i>Recommended Change:</i>	<i>12 months</i>
LC195	Optipure L1 Purification Cartridge
<i>Recommended Change:</i>	<i>3-12 months</i>
LC202	Vent Filter
<i>Recommended Change:</i>	<i>6 months</i>
CT3-50	Pack of 50 Effervescent Sanitisation Tablets
CT3-300	6 Packs of 50 Effervescent Sanitisation Tablets
<i>Typical usage:</i>	<i>10 tablets per month</i>

## MEDICA Biopure 60/120 Accessories

LA682	Hygienic sample point (15mm & 22mm)
LA701	E-Cartridge Upgrade Kit complete with LC181
LA703	60LPH to 120LPH Upgrade Kit
LA704	Castors



**This guide is not for use in USA**

# BIOPURE

Part Number	Description
MB200XXM1-230	BIOPURE 200
MB300XXM1-230	BIOPURE 300
MB600XXM1-230	BIOPURE 600

Unit is supplied as standard with:

Fitted / Included: 1 x LA682 and sampling accessory bag, Appropriate type and qty RO module 1 x LC157, 1 x LC158, Operator manual, PASSkey set, 1 x CFYAZD9960, 2x Filter spanners

## MEDICA Biopure 200/300/600 Consumables

LC157	Bacteriological reservoir filter
<i>Recommended Change:</i>	<i>6 months</i>
LC161	Reverse Osmosis Module
LC176	Reverse Osmosis Cartridge 300lph
<i>Recommended Change:</i>	<i>Typically 2-3 years</i>
LC158	UV Lamp
<i>Recommended Change:</i>	<i>12 months</i>
LC172	Big Blue 20" Carbon Filter
<i>Recommended Change:</i>	<i>6 months</i>
CFYAZD9960	Ultra-Filter
<i>Recommended Change:</i>	<i>12 months</i>

## MEDICA Biopure 200/300/600 Accessories

LA665	Remote Display
LA667	Replacement Pass Key (Blue Operator Key)
LA668	Replacement Pass Key (Green Sanitisation Key)
LA681	Filter Bracket Assembly
LA682	Hygienic sample point (15mm & 22mm)
LA683	Filter Housing - Big Blue
LA688	Temperature Control Kit
LA679	High recovery (75%) Reverse Osmosis kit
LA692	High Recovery Kit

## Aquadem Cylinder (example)

PRDISI0067	E 1000 Cylinder, with 94 litres of MB20 resins, fitted with High Flow Head
VMGUVA200684	High flow top cap with 3/4" Male Connectors
FPPIP203045	Tubing set for use with High Flow Top Cap. 3/4" female Connectors



**This guide is not for use in USA**

# CENTRA R 60/120

Part Number	Description
CN060RDM1-115	CENTRA R 60 115V
CN060RDM1-230	CENTRA R 60 230V
CN120RDM1-115	CENTRA R 120 115V
CN120RDM1-230	CENTRA R 120 230V

Unit is supplied as standard with:

Fitted / Included: 1 x LC118, 2 x LC136M2, Appropriate RO type and Quantity for unit, Quick Start reference and operator manuals, user recognition key set, LA724 bypass block

## CENTRA R 60 /120 Starter Kits

LA721	Starter Kit 1
Contains	1 x LC175, LA720 installation kit
LA722	Starter Kit 2
Contains	1 x LC177, LA720 installation kit

## CENTRA R 60 /120 Consumables

LC175	Protek L1 Pretreatment Cartridge
Recommended Change:	6 months
LC177	Protek L2 Pretreatment Cartridge
Recommended Change:	12 months
LC180	60L/Hr Reverse Osmosis Cartridge
Recommended Change:	Typically 3 Years
LC181	E - Cartridge
Recommended Change:	3 Years
LC118	UV Lamp
Recommended Change:	12 months
LC187	Purification Pack
Recommended Change:	3-12 months
LC136M2	Composite Air Vent Filter
Recommended Change	(2 of): 12 months
LC204	Microfilter
Recommended Change:	12 months
CT3	Pack of 12 Effervescent Sanitisation Tablets
CT3-50	Pack of 50 Effervescent Sanitisation Tablets
CT3-300	6 Packs of 50 Effervescent Sanitisation Tablets
Typical usage:	6 tablets per month

## CENTRA R 60 /120 Accessories:

LA701	E-Cartridge Upgrade Kit complete with LC181
LA703	60LPH to 120LPH Upgrade Kit
LA704	Castors
LA631	40 Litre Tank complete with Level Control Assembly 230V 50Hz
LA633	40 Litre Tank complete with Level Control Assembly 115V 60Hz
LA643	Remote Dispense Station (Including Dispense Gun)



**This guide is not for use in USA**

# CENTRA R200

Part Number	Description
CN200RDM1-115	CENTRA R 200 (115V 60Hz) 200 LPH Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration
CN200RDM1-230	CENTRA R 200 (230V 50Hz) 200 LPH Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration
CNHFR2M1-230	CENTRA R 200 High Recirculation (230V 50Hz) 200 LPH Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration
CNHFVR2M1-230	CENTRA R 200 Variable Speed (230V 50Hz) 200 LPH Reverse Osmosis/Central Distribution System complete with 350 Litre Tank, UV & Filtration

Unit is supplied as standard with:

Fitted / Included: 1 x LC156/LC173 x LC159, 1 x LC161, 1 x LC158, 1 x LC160 (non-high flow models only), Operator manual, PASSkey set, 2 x Filter spanners and 2 x sanitisation tubes

## Consumables

LC159	Pre-treatment (pump protection) Filter
Recommended Change:	6 months
LC156	Composite Vent Filter
Recommended Change:	6 months
LC173	Composite Vent Filter
Recommended Change:	6 months
LC158	UV Lamp
Recommended Change:	12 months
LC160	0.2 micron filter
Recommended Change:	6 months

## Accessories

LA665	Remote Display
LA667	Replacement Pass Key (Blue Operator Key)
LA668	Replacement Pass Key (Green Sanitisation Key)
LA679	High recovery (75%) Reverse Osmosis kit

## Aquadem Cylinder (example)

PRDISI0075	E 1000 Cylinder, with 94 litres of HYPEX resins, fitted with High Flow Head
VMGUVA200684	High flow top cap with 3/4" Male Connectors
FPPIPV203045	Tubing set for use with High Flow Top Cap. 3/4" female Connectors



**This guide is not for use in USA**

