

The right temperature worldwide

**LAUDA**

how  
**PRO**  
are you?



Discover the new benchmark for no-compromise temperature control:  
**LAUDA PRO bath and circulation thermostats**



# The new LAUDA PRO: compact, efficient, powerful. Optimized to your application.

## PRO or PRO? Your choice.

Real professionals do not make compromises – least of all with temperature control. We developed the new LAUDA PRO thermostats for you.

Powerful bath thermostats for precise internal temperature control from -100 up to 250 °C. And particularly efficient circulation thermostats for dynamic heating and cooling of external applications from -45 up to 250 °C.

Professional laboratory technicians get what they really need with the new LAUDA PRO: no-compromise reliable temperature control in the research environment, application technology and production. **How PRO are you?**



**100 % flexibility**  
**for internal temperature**  
**control.**



**Application examples**

Functional check of electronic components  
Material testing of metal and plastic samples  
Sensor calibration  
Temperature control of measuring cells

**100 % efficiency**  
**for external applications.**



**Application examples**

Thermal control of chemical and biological reactions  
Precise temperature control for crystallization processes  
Temperature control of measuring and testing rigs, heat exchangers and distillation plant



- + Rapid temperature change
- + Lower operating costs thanks to low energy and material consumption
- + Control possible from up to 50 meters away
- + Waste heat reduction as a result of hybrid cooling
- + USB and Ethernet interfaces are standard and other interface modules are available

Your choice:

## BASE control unit.

### Flexible and user-friendly.

Control your LAUDA PRO just by pushing a button and maintain an overview on the brilliant OLED displays.

#### Your advantages at a glance

- + Bright high contrast OLED display
- + Convenient control using softkeys and cursor keys
- + Clear menu navigation
- + Programmer with 1 program and 20 segments
- + Removable
- + Control unit for the bath thermostats can rotate through 360°
- + Temperature display in °C or °F
- + Cable extension up to 50 meters



Your choice:

## COMMAND control unit.

### Convenient multi-touch screen.

Even more convenience is offered by the COMMAND control unit, with its large color multi-touch screen and versatile functions.

#### Your advantages at a glance

- + Large color touch screen (5.7")
- + Multi-touch operation using freely configurable tiles – can also be operated wearing gloves
- + Zoomable graphic view
- + Programmer with 100 programs each having 50 segments
- + Weekly timer and count-down function
- + Ramp function
- + Data logging, export to USB stick
- + Individual user rights for up to 20 users
- + Removable
- + Control unit for the bath thermostats can rotate through 360°
- + Temperature display in °C or °F
- + Cable extension up to 50 meters



Reliable operating convenience - with and without gloves.





**Our product range:**

Thermostats · Circulation Chillers · Water Baths

Process Cooling Units · Heat Transfer Systems · Secondary Circuit Units

**For professionals who want to know more:**  
**how-pro-are-you.com**



Further information and a practical product finder  
await you at the website:

**[www.how-pro-are-you.com](http://www.how-pro-are-you.com)**



Your dealer

LAUDA DR. R. WOBSE GMBH & CO. KG  
Pfarrstrasse 41/43 · 97922 Lauda-Königshofen · Germany  
Phone.: +49 (0)9343 503-0 · Fax: +49 (0)9343 503-222  
Email: [info@lauda.de](mailto:info@lauda.de) · Internet: [www.lauda.de](http://www.lauda.de)