

## MTPS Series with Flex Option

# MULTI-CONFIGURABLE PLANT GROWTH ROOM

Adaptable to a Range of Research Programs



## With the Flex option, adapt to a range of plant heights, lighting intensities and airflow to suit research requirements

#### **TEMPERATURE**

 +10°C to +35°C Lights OFF +15°C to +35°C Lights ON

#### **REFRIGERATION**

• Self-contained water-cooled condensing unit

#### **MAINTENANCE**

• Easy access to refrigeration and humidification equipment through partially open top and removable upper front fascia panel

#### **HUMIDITY**

- Optional additive humidity
- Up to 70%RH\*
- Control accuracy: ±3%\*

#### **ACCESS**

- Clear opening door with inside safety release latch and cam-type selfclosing hinges.
- Door has positive closer device, thermal plastic gasket with magnetic core, and door stop

### **GROWTH HEIGHTS**

- One tier: 731/2" (1865 mm)
- Two tiers: 36" (914 mm)
- Three tiers: 20" (510 mm)


#### LIGHTING

- Fluorescent lighting up to:
  - » One tier: 1000 µmol/m<sup>2</sup>/s
  - » Two tier: 700 µmol/m²/s per tier
  - » Three tier: 500 µmol/m²/s per tier




#### **AIRFLOW**


Adjusts with a simple flip of a switch



**SINGLE TIER** UPWARD or HORIZONTAL AIRFLOW (Upward Shown)



**TWO TIER** HORIZONTAL AIRFLOW



THREE TIER HORIZONTAL AIRFLOW

July 2017, Rev 02 | MK0051



<sup>\*</sup>Consult factory for details