

Thermo Scientific™  
Finnpipette™ Novus™  
Electronic Pipettes and  
Finntip™ Tips

- 10 µL
- 20 µL
- 50 µL
- 100 µL /200 µL
- 300 µL
- 1000 µL
- 1200 µL\*
- 5 mL
- 10 mL

\*Finntip Flex tips only

| Color Key/Tip Packaging/Volume Range                             |                   |                 |                          |
|------------------------------------------------------------------|-------------------|-----------------|--------------------------|
| 10 µL MICRO PIPETTES                                             | Finnpipette Novus | Increment       | Functional Volume Range¹ |
| • 1-10 µL                                                        | 46200000          | 0.01 µL         | 0.5-10 µL                |
| • 8-ch 1-10 µL                                                   | 46300000          | 0.01 µL         | 0.5-10 µL                |
| • 12-ch 1-10 µL                                                  | 46300100          | 0.01 µL         | 0.5-10 µL                |
| 10 µL FINNTIPS                                                   | Finntip 10        | Finntip Flex 10 | Tip Volume Range         |
| • Racked, 96 tips/rack, 10 racks/pack                            | 9400300           | 94060100        | 0.2 - 10 µL              |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                   | 9400303           | 94060103        | 0.2 - 10 µL              |
| • Reload Starter Kit, 4 x 192 tips/rack + 10 x 192 tips/insert   | 9400326           | 94060116*       | 0.2 - 10 µL              |
| • Reload Tower, 192 tips/insert, 20 inserts/pack                 | 9400327           | 94060117**      | 0.2 - 10 µL              |
| • Bulk, 1000 tips/bag                                            | 9400310           | 94060120        | 0.2 - 10 µL              |
| • Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94052000          | 94056980        | 0.2 - 10 µL              |

\*Finntip Flex 10 Reload Starter Kit includes 1 x 96 tips/rack + 20 x 96 tips/insert

\*\*Finntip Flex Reload Kit includes 96 tips/insert, 20 inserts/pack

| 50 µL MICRO PIPETTES                                                     | Finnpipette Novus | Increment        | Functional Volume Range¹ |
|--------------------------------------------------------------------------|-------------------|------------------|--------------------------|
| • 5-50 µL                                                                | 46200200          | 0.1 µL           | 2.5-50 µL                |
| • 16-ch 5-50 µL                                                          | 46300700          | 0.1 µL           | 2.5-50 µL                |
| 50 µL FINNTIPS MICRO, 384-FORMAT FINNTIPS                                | Finntip 50        | Tip Volume Range |                          |
| • Racked, 384 tips/rack, 10 racks/pack                                   | 9400370           | 0.2 - 50 µL      |                          |
| • Racked, Sterile, 384 tips/rack, 10 racks/pack                          | 9400373           | 0.2 - 50 µL      |                          |
| • Bulk, 1000 tips/bag                                                    | 9400360           | 0.2 - 50 µL      |                          |
| • Racked, <b>Filtered</b> , Sterile, Micro, 384 tips/rack, 10 racks/pack | 94052060          | 0.2 - 50 µL      |                          |

| 10/50/100 µL PIPETTES | Finnpipette Novus | Increment | Functional Volume Range¹ |
|-----------------------|-------------------|-----------|--------------------------|
| • 1-10 µL             | 46200100          | 0.01 µL   | 0.5-10 µL                |
| • 5-50 µL             | 46200300          | 0.1 µL    | 2.5-50 µL                |
| • 10-100 µL           | 46200400          | 0.1 µL    | 5-100 µL                 |
| • 8-ch 5-50 µL        | 46300200          | 0.1 µL    | 2.5-50 µL                |
| • 12-ch 5-50 µL       | 46300300          | 0.1 µL    | 2.5-50 µL                |

| 200/250 µL FINNTIPS                                          | Finntip 250 | Tip Volume Range | Finntip Flex 200 | Tip Volume Range |
|--------------------------------------------------------------|-------------|------------------|------------------|------------------|
| • Racked, 96 tips/rack, 10 racks/pack                        | 9400260     | 0.5 - 250 µL     | 94060310         | 1 - 200 µL       |
| • Racked, sterile, 96 tips/rack, 10 racks/pack               | 9400263     | 0.5 - 250 µL     | 94060313         | 1 - 200 µL       |
| • Reload Starter Kit, 4 x 96 tips/rack + 10 x 96 tips/insert | 9400266     | 0.5 - 250 µL     | 94060316*        | 1 - 200 µL       |
| • Reload Tower, 96 tips/insert, 20 inserts/pack              | 9400267     | 0.5 - 250 µL     | 94060317         | 1 - 200 µL       |
| • Bulk, 500 tips/box                                         | 9400250     | 0.5 - 250 µL     | —                | —                |
| • Bulk, 1000 tips/bag                                        | 9400230     | 0.5 - 250 µL     | 94060320         | 1 - 200 µL       |

\*Finntip Flex 200 Reload Starter Kit includes 1 x 96 tips/rack + 20 x 96 tips/insert

| 200/250 µL FINNTIP Extended Length and Wide  | Finntip 200 Ext <sup>2</sup> Tip Volume Range |            | Finntip 250 Wide <sup>3</sup> | Tip Volume Range |
|----------------------------------------------|-----------------------------------------------|------------|-------------------------------|------------------|
| Racked, 96 tips/rack, 10 racks/pack          | 9400130                                       | 5 - 200 µL | 9405120                       | 10 - 250 µL      |
| Racked, sterile, 96 tips/rack, 10 racks/pack | 9400133                                       | 5 - 200 µL | 9405123                       | 10 - 250 µL      |
| Bulk, 500 tips/box                           | 9400100**                                     | 5 - 200 µL | —                             | —                |
| Bulk, 1000 tips/bag                          | —                                             | —          | 9405020                       | 10 - 250 µL      |

\*\*Finntip 200 Ext Bulk pack includes 400 tips/box

¹The Functional volume range indicates the volume range that the pipette can cover in the Stepper, Sequential Stepper and Sequential Aspirate functions.

²Color coded orange.

³Color coded white

| Color Key/Tip Packaging/Volume Range                                                |          |                  |
|-------------------------------------------------------------------------------------|----------|------------------|
| 10/100/200 µL FILTERED FINNTIPS                                                     | Finntip  | Tip Volume Range |
| • Finntip 10 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack       | 94052100 | 0.5 - 10 µL      |
| • Finntip 100 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack      | 94052200 | 0.5 - 100 µL     |
| • Finntip 100 Ext / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack  | 94052310 | 5 - 100 µL       |
| • Finntip Flex 100 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94056520 | 1 - 100 µL       |
| • Finntip Flex 200 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94056380 | 1 - 200 µL       |

The universal-fit yellow color coded filtered tips are marked differently from other Finntips. The filter is in a different location based on tip size. We recommend matching the filtered tip size with the pipette volume range in addition to color code.

| 300 µL PIPETTES                                                 | Finnpipette Novus | Increment        | Functional Volume Range¹ |
|-----------------------------------------------------------------|-------------------|------------------|--------------------------|
| • 30-300 µL                                                     | 46200500          | 1 µL             | 15-300 µL                |
| • 8-ch 30-300 µL                                                | 46300400          | 1 µL             | 15-300 µL                |
| • 12-ch 30-300 µL                                               | 46300500          | 1 µL             | 15-300 µL                |
| 300 µL FINNTIPS                                                 | Finntip 300       | Finntip Flex 300 | Tip Volume Range         |
| • Racked, 96 tips/rack, 10 racks/pack                           | 9401250           | 94060510         | 5 - 300 µL               |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                  | 9401253           | 94060513         | 5 - 300 µL               |
| • Reload Starter Kit 1 x 96 tips/rack +20 x 96 tips/insert      | —                 | 94060516         | 5 - 300 µL               |
| • Reload Tower, 96 tips/insert, 20 inserts/pack                 | 9401255*          | 94060517         | 5 - 300 µL               |
| • Bulk, 1000 tips/bag                                           | 9401260           | 94060520         | 5 - 300 µL               |
| • Racked, <b>Filtered</b> , sterile, 96 tip/rack, 10 racks/pack | 94052350          | 94056580         | 5 - 300 µL               |

\* Finntip 300 Reload package includes 10 units of 2 x 96 tips/insert

| 1000 µL PIPETTES                                                 | Finnpipette Novus |                   | Increment        | Functional Volume Range¹ |                  |
|------------------------------------------------------------------|-------------------|-------------------|------------------|--------------------------|------------------|
| • 100-1000 µL                                                    | 46200600          |                   | 1 µL             | 50-1000 µL               |                  |
| 1000 µL FINNTIPS                                                 | Finntip 1000      | Finntip Flex 1000 | Finntip 1000 Ext | Finntip Wide 1000³       | Tip Volume Range |
| • Racked, 96 tips/rack, 10 racks/pack                            | 9401110           | 94060710          | 9401420          | 9405160                  | 100-1000 µL      |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                   | 9401113           | 94060713          | 9401423          | 9405163                  | 100-1000 µL      |
| • Reload Starter Kit 1 x 96 tips/rack + 16 x 96 tips/insert      | —                 | 94060716          | —                | —                        | 100-1000 µL      |
| • Reload Tower, 96 tips/insert, 16 inserts/pack                  | 9401115*          | 94060717          | —                | —                        | 100-1000 µL      |
| • Bulk, 200 tips/box                                             | 9401070           | —                 | —                | —                        | 100-1000 µL      |
| Bulk, 400 tips/bag                                               | —                 | —                 | —                | 9405050                  | 100-1000 µL      |
| • Bulk, 1000 tips/bag                                            | 9401030           | 94060720          | 9401410          | —                        | 100-1000 µL      |
| • Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94052410          | 94056710          | 94052430         | —                        | 100-1000 µL      |

\*Finntip 1000 Reload package includes 5 units of 2 x 96 tips/insert

| 1200 µL PIPETTES                                                | Finnpipette Novus | Increment        | Functional Volume Range¹ |
|-----------------------------------------------------------------|-------------------|------------------|--------------------------|
| • 8-ch 100-1200 µL                                              | 46300800          | 1 µL             | 50-1200 µL               |
| 1200 µL FINNTIP FLEX TIPS                                       | Finntip Flex 1200 | Tip Volume Range |                          |
| • Racked, 96 tips/rack, 10 racks/pack                           | 94060810          | 50 - 1200 µL     |                          |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                  | 94060813          | 50 - 1200 µL     |                          |
| • Reload Starter Kit 1 rack 96 tips/rack + 16 x 96 tips/inserts | 94060816          | 50 - 1200 µL     |                          |
| • Reload Kit, 16 x 96 tips/insert                               | 94060817          | 50 - 1200 µL     |                          |
| • Bulk, 1000 tips/bag                                           | 94060820          | 50 - 1200 µL     |                          |
| • Racked, <b>Filtered</b> , sterile, 96tips/rack, 10 racks/pack | 94056810          | 50 - 1200 µL     |                          |

| 5 mL PIPETTES                                                   | Finnpipette Novus | Increment        | Functional Volume Range¹ |
|-----------------------------------------------------------------|-------------------|------------------|--------------------------|
| • 0.5-5 mL                                                      | 46200700          | 0.01 mL          | 0.25 - 5 mL              |
| 5 mL FINNTIPS                                                   | Finntip 5 mL      | Tip Volume Range |                          |
| • Racked, 54 tips/rack, 5 racks/pack                            | 9402070           | 0.5 - 5 mL       |                          |
| • Racked, sterile, 54 tips/rack, 5 racks/pack                   | 9402073           | 0.5 - 5 mL       |                          |
| • Bulk, 75 tips/box                                             | 9402050           | 0.5 - 5 mL       |                          |
| • Bulk, 5 x 100 tips/bag                                        | 9402030           | 0.5 - 5 mL       |                          |
| • Bulk box, 30 x 100 tips/bag                                   | 9402010           | 0.5 - 5 mL       |                          |
| • Racked, <b>Filtered</b> , sterile, 54 tips/rack, 5 racks/pack | 94052550          | 0.5 - 5 mL       |                          |

| 10 mL PIPETTES                                                  | Finnpipette Novus | Increment              | Functional Volume Range¹ |
|-----------------------------------------------------------------|-------------------|------------------------|--------------------------|
| • 1-10 mL                                                       | 46200800          | 0.01 mL                | 0.5-10 mL                |
| 10 mL FINNTIPS                                                  | Finntip 10 mL     | Finntip Flex 10 mL Ext | Tip Volume Range         |
| • Racked, 24 tips/rack, 5 racks/pack                            | 9402160           | —                      | 1 - 10 mL                |
| • Racked, sterile, 24 tips/rack, 5 racks/pack                   | 9402163           | —                      | 1 - 10 mL                |
| • Bulk, 40 tips/box                                             | 9402150           | —                      | 1 - 10 mL                |
| • Bulk, 100 tips/bag                                            | 9402151           | 94060970               | 1 - 10 mL                |
| • Racked, <b>Filtered</b> , sterile, 24 tips/rack, 5 racks/pack | 94052600          | 94056970               | 1 - 10 mL                |

Thermo Scientific™ Finnpipette Novus Good Laboratory Pipetting (GLP) KITS

| Kit Number | Content                                                                                                                                                       | Item Number |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| GLP Kit 1  | • Single channel 10-100 µL and 100-1000 µL pipettes<br>• Two Novus Stands<br>• Finntip Flex 200 and Flex 1000 racks                                           | 4700470     |
| GLP Kit 2  | • Single channel 100-1000 and 0.5-5 mL pipettes<br>• Two Novus Stands<br>• Finntip Flex 1000 and Finntip 5 mL racks                                           | 4700490     |
| GLP Kit 3  | • Single channel 1-10 µL micro, 10-100 µL and 100-1000 µL pipettes<br>• Three Novus Stands<br>• Finntip Flex 10, Flex 200 and Flex 1000 racks                 | 4701120     |
| GLP Kit 4  | • Single channel 10-100 µL, 100-1000 µL and 1-10 mL pipettes<br>• Three Novus Stands<br>• Finntip Flex 200, Flex 1000, Finntip 10 mL and Flex 10 mL Ext racks | 4701130     |



Stands for Manual and Electronic Finnpipettes

| Description                                                               | Item Number |
|---------------------------------------------------------------------------|-------------|
| F-Series Stand, white, linear, for six F1 or F2 pipettes                  | 9420400     |
| Multichannel Pipette Stand, for one F1, F2 or Novus pipette               | 9420390     |
| Pipette Stand, grey, linear, for six F1 or F2 pipettes                    | 9420290     |
| Pipette Mini-stand, linear, for three F1 or F2 pipettes                   | 9420320     |
| Carousel Stand, for six F1 or F2, or three Novus pipettes, includes cover | 9420340     |
| Novus stand, for one pipette                                              | 9420360     |



F-Series Stand



Multichannel Pipette Stand



Carousel Stand



Novus Stand

## Thermo Scientific Finnpipette F1, F2 and Novus Pipette Systems

Selection Guide

Find out more at [thermofisher.com/finnpipette](https://thermofisher.com/finnpipette)

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Please consult your local sales representative for details.

ThermoFisher  
SCIENTIFIC



ThermoFisher  
SCIENTIFIC



# Thermo Scientific Finnpipette F1 and F2 Manual Pipettes and Finntip Tips

- 10 µL
- 20 µL
- 50 µL
- 100 µL/200 µL
- 300 µL
- 1000 µL
- 1200 µL\*
- 5 mL
- 10 mL

\*Finntip Flex tips only

| Color Key/Tip Packaging/Volume Range                             |                |                 |                  |
|------------------------------------------------------------------|----------------|-----------------|------------------|
| 2/10 µL MICRO PIPETTES                                           | Finnpipette F1 | Finnpipette F2  | Increment        |
| • 0.2-2 µL                                                       | 4641010N       | 4642010         | 0.002 µL         |
| • 0.5-5 µL                                                       | 4641020N       | 4642020         | 0.01 µL          |
| • 1-10 µL                                                        | 4641030N       | 4642030         | 0.02 µL          |
| • 8-ch 1-10 µL                                                   | 4661000N       | 4662000         | 0.02 µL          |
| • 12-ch 1-10 µL                                                  | 4661040N       | 4662040         | 0.02 µL          |
| 10 µL FINNTIPS                                                   | Finntip 10     | Finntip Flex 10 | Tip Volume Range |
| • Racked, 96 tips/rack, 10 racks/pack                            | 9400300        | 94060100        | 0.2 - 10 µL      |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                   | 9400303        | 94060103        | 0.2 - 10 µL      |
| • Reload Starter Kit, 4 x 192 tips/rack + 10 x 192 tips/insert   | 9400326        | 94060116*       | 0.2 - 10 µL      |
| • Reload Tower, 192 tips/insert, 20 inserts/pack                 | 9400327        | 94060117**      | 0.2 - 10 µL      |
| • Bulk, 1000 tips/bag                                            | 9400310        | 94060120        | 0.2 - 10 µL      |
| • <b>Filtered</b> , Racked, sterile, 96 tips/rack, 10 racks/pack | 94052000       | 94056980        | 0.2 - 10 µL      |

\*Finntip Flex 10 Reload Starter Kit includes 1 x 96 tips/rack + 20 x 96 tips/insert

\*\*Finntip Flex Reload Tower includes 96 tips/insert, 20 inserts/pack

| 10 µL MICRO, 384-FORMAT PIPETTES                                  | Finnpipette F1 | Finnpipette F2 | Increment        |
|-------------------------------------------------------------------|----------------|----------------|------------------|
| • 16-ch 1-10 µL                                                   | 4661080N       | 4662080        | 0.02 µL          |
| 20 µL FINNTIP MICRO, 384-FORMAT TIPS                              |                | Finntip 20     | Tip Volume Range |
| • Racked, 384 tips/rack, 10 racks/pack                            |                | 9400610        | 0.2-20 µL        |
| • Racked, Sterile, 384 tips/rack, 10 racks/pack                   |                | 9400613        | 0.2-20 µL        |
| • Bulk, 1000 tips/bag                                             |                | 9400620        | 0.2-20 µL        |
| • <b>Filtered</b> , Racked, Sterile, 384 tips/rack, 10 racks/pack |                | 94052020       | 0.2-20 µL        |

| 20/50 µL MICRO & 384-FORMAT PIPETTES                              | Finnpipette F1 | Finnpipette F2 | Increment        |
|-------------------------------------------------------------------|----------------|----------------|------------------|
| • 2-20 µL                                                         | 4641050N       | 4642050        | 0.02 µL          |
| • 5-50 µL                                                         | 4641130N       | 4642120        | 0.1 µL           |
| • 16-ch 5-50 µL                                                   | 4661090N       | 4662090        | 0.1 µL           |
| 50 µL FINNTIP MICRO, 384-FORMAT TIPS                              |                | Finntip 50     | Tip Volume Range |
| • Racked, 384 tips/rack, 10 racks/pack                            |                | 9400370        | 0.2-50 µL        |
| • Racked, Sterile, 384 tips/rack, 10 racks/pack                   |                | 9400373        | 0.2-50 µL        |
| • Bulk, 1000 tips/bag                                             |                | 9400360        | 0.2-50 µL        |
| • <b>Filtered</b> , Racked, Sterile, 384 tips/rack, 10 racks/pack |                | 94052060       | 0.2-50 µL        |

| 10/200 µL PIPETTES | Finnpipette F1 | Finnpipette F2 | Increment |
|--------------------|----------------|----------------|-----------|
| • 1-10 µL          | 4641040N       | 4642040        | 0.02 µL   |
| • 2-20 µL          | 4641060N       | 4642060        | 0.02 µL   |
| • 5-50 µL          | 4641140N       | 4642130        | 0.1 µL    |
| • 10-100 µL        | 4641070N       | 4642070        | 0.2 µL    |
| • 20-200 µL        | 4641080N       | 4642080        | 0.2 µL    |
| • 8-ch 5-50 µL     | 4661010N       | 4662010        | 0.1 µL    |
| • 8-ch 10-100 µL   | 4661020N       | 4662020        | 0.2 µL    |
| • 12-ch 5-50 µL    | 4661050N       | 4662050        | 0.1 µL    |
| • 12-ch 10-100 µL  | 4661060N       | 4662060        | 0.2 µL    |

| Color Key/Tip Packaging/Volume Range                         |                                               |                                                |             |
|--------------------------------------------------------------|-----------------------------------------------|------------------------------------------------|-------------|
| 200/250 µL FINNTIPS                                          | Finntip 250 Tip Volume Range                  | Finntip Flex 200 Tip Volume Range              |             |
| • Racked, 96 tips/rack, 10 racks/pack                        | 9400260                                       | 94060310                                       | 1 - 200 µL  |
| • Racked, sterile, 96 tips/rack, 10 racks/pack               | 9400263                                       | 94060313                                       | 1 - 200 µL  |
| • Reload Starter Kit, 4 x 96 tips/rack + 10 x 96 tips/insert | 9400266                                       | 94060316*                                      | 1 - 200 µL  |
| • Reload Tower, 96 tips/insert, 20 inserts/pack              | 9400267                                       | 94060317                                       | 1 - 200 µL  |
| • Bulk, 500 tips/box                                         | 9400250                                       | –                                              | –           |
| • Bulk, 1000 tips/bag                                        | 9400230                                       | 94060320                                       | 1 - 200 µL  |
| 200/250 µL FINNTIP EXTENDED LENGTH AND WIDE TIPS             | Finntip 200 Ext <sup>1</sup> Tip Volume Range | Finntip 250 Wide Tip Volume Range <sup>2</sup> |             |
| Racked, 96 tips/rack, 10 racks/pack                          | 9400130                                       | 9405120                                        | 10 - 250 µL |
| Racked, sterile, 96 tips/rack, 10 racks/pack                 | 9400133                                       | 9405123                                        | 10 - 250 µL |
| Bulk, 500 tips/box                                           | 9400100**                                     | –                                              | 10 - 250 µL |
| Bulk, 1000 tips/bag                                          | –                                             | 9405020                                        | 10 - 250 µL |

\*Finntip Flex 200 Reload Starter Kit includes 1 x 96 tips/rack + 20 x 96 tips/insert

\*\*Finntip 200 Ext Bulk pack includes 400 tips/box

| Filtered FINNTIPS 10-200 2/10 µL                                                    | Finntip               | Tip Volume Range |
|-------------------------------------------------------------------------------------|-----------------------|------------------|
| • Finntip 10 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack       | 94052100              | 0.5 - 10 µL      |
| • Finntip 20 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack       | 94052150              | 0.5 - 20 µL      |
| • Finntip 30 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack       | 94052160              | 0.5 - 30 µL      |
| • Finntip 100 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack      | 94052200              | 0.5 - 100 µL     |
| • Finntip 100 Ext / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack  | 94052310 <sup>1</sup> | 5 - 100 µL       |
| • Finntip 200 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack      | 94052300              | 0.5 - 200 µL     |
| • Finntip 200 Ext / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack  | 94052320 <sup>1</sup> | 5 - 200 µL       |
| • Finntip Flex 30 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack  | 94056510              | 1 - 30 µL        |
| • Finntip Flex 100 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94056520              | 1 - 100 µL       |
| • Finntip Flex 200 / Racked, <b>Filtered</b> , sterile, 96 tips/rack, 10 racks/pack | 94056380              | 1 - 200 uL       |

The universal-fit yellow color coded filtered tips are marked differently from other Finntips. The filter is at different location based on tip size. Match the filtered tip size with the pipette volume range in addition to color code.

| 300 µL PIPETTES                                                 | Finnpipette F1 | Finnpipette F2   | Increment        |
|-----------------------------------------------------------------|----------------|------------------|------------------|
| • 30-300 µL                                                     | 4641090N       | –                | 1 µL             |
| • 8-ch 30-300 µL                                                | 4661030N       | 4662030          | 1 µL             |
| • 12-ch 30-300 µL                                               | 4661070N       | 4662070          | 1 µL             |
| 300 µL FINNTIPS                                                 | Finntip 300    | Finntip Flex 300 | Tip Volume Range |
| • Racked, 96 tips/rack, 10 racks/pack                           | 9401250        | 94060510         | 5 - 300 µL       |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                  | 9401253        | 94060513         | 5 - 300 µL       |
| • Reload Starter Kit 1 x 96 tips/rack + 20 x 96 tips/insert     | –              | 94060516         | 5 - 300 µL       |
| • Reload Tower, 96 tips/insert, 20 inserts/pack                 | 9401255*       | 94060517         | 5 - 300 µL       |
| • Bulk, 1000 tips/bag                                           | 9401260        | 94060520         | 5 - 300 µL       |
| • <b>Filtered</b> , Racked, sterile, 96 tip/rack, 10 racks/pack | 94052350       | 94056580         | 5 - 300 µL       |

\* Finntip 300 Reload package includes 10 units of 2 x 96 tips/insert

| 1000 µL PIPETTES                                                 | Finnpipette F1 | Finnpipette F2    | Increment           |                                |
|------------------------------------------------------------------|----------------|-------------------|---------------------|--------------------------------|
| • 100-1000 µL                                                    | 4641100N       | 4642090           | 1 µL                |                                |
| 1000 µL FINNTIPS                                                 | Finntip 1000   | Finntip Flex 1000 | Finntip 1000 µL Ext | Finntip Wide 1000 <sup>2</sup> |
| • Racked, 96 tips/rack, 10 racks/pack                            | 9401110        | 94060710          | 9401420             | 9405160                        |
| • Racked, sterile, 96 tips/rack, 10 racks/pack                   | 9401113        | 94060713          | 9401423             | 9405163                        |
| • Reload Starter Kit 1 x 96 tips/rack + 16 x 96 tips/insert      | —              | 94060716          | —                   | —                              |
| • Reload Tower, 96 tips/insert, 16 inserts/pack                  | 9401115*       | 94060717          | —                   | —                              |
| • Bulk, 200 tips/box                                             | 9401070        | —                 | —                   | —                              |
| • Bulk, 400 tips/bag                                             | —              | —                 | —                   | 9405050                        |
| • Bulk, 1000 tips/bag                                            | 9401030        | 94060720          | 9401410             | —                              |
| • <b>Filtered</b> , racked, sterile, 96 tips/rack, 10 racks/pack | 94052410       | 94056710          | 94052430            | —                              |

\* Finntip 1000 Reload package includes 5 units of 2 x 96 tips/insert

<sup>1</sup> Color coded orange

<sup>2</sup> Color coded white

| Color Key/Tip Packaging/Volume Range                            |                |                |                  |
|-----------------------------------------------------------------|----------------|----------------|------------------|
| 5 mL PIPETTES                                                   | Finnpipette F1 | Finnpipette F2 | Increment        |
| • 0.5-5 mL                                                      | 4641110N       | 4642100        | 0.01 mL          |
| 5 mL FINNTIPS                                                   |                | Finntip 5 mL   | Tip Volume Range |
| • Racked, 54 tips/rack, 5 racks/pack                            |                | 9402070        | 0.5 - 5 mL       |
| • Racked, sterile, 54 tips/rack, 5 racks/pack                   |                | 9402073        | 0.5 - 5 mL       |
| • Bulk, 5 x 100 tips/bag                                        |                | 9402030        | 0.5 - 5 mL       |
| • Bulk box, 30 x 100 tips/bag                                   |                | 9402 010       | 0.5 - 5 mL       |
| • Bulk, 75 tips/box                                             |                | 9402050        | 0.5 - 5 mL       |
| • <b>Filtered</b> , racked, sterile, 54 tips/rack, 5 racks/pack |                | 94052550       | 0.5 - 5 mL       |

| 10 mL PIPETTES                                                   | Finnpipette F1 | Finnpipette F2        | Increment        |
|------------------------------------------------------------------|----------------|-----------------------|------------------|
| • 1-10 mL                                                        | 4641120N       | 4642110               | 0.02 mL          |
| 10 mL FINNTIPS                                                   | Finntip 10 mL  | Finntip Flex 10mL Ext | Tip Volume Range |
| • Racked, 24 tips/rack, 5 racks/pack                             | 9402160        | –                     | 1 - 10 mL        |
| • Racked, sterile, 24 tips/rack, 5 racks/pack                    | 9402163        | –                     | 1 - 10 mL        |
| • Bulk, 100 tips/bag                                             | 9402151        | 94060970              | 1 - 10 mL        |
| • Bulk, 40 tips/box                                              | 9402150        | –                     | 1 - 10 mL        |
| • <b>Filtered</b> , racked, sterile, 24 tips/rack, 5 raacks/pack | 94052600       | 94056970              | 1 - 10 mL        |

## Finnpipette F1 and F2 Fixed Volume Pipettes

| Volume   | Finnpipette F1 | Finnpipette F2 |
|----------|----------------|----------------|
| 1 µl     | •              | 4652000        |
| 5 µl     | –              | 4652010        |
| 10 µl    | 4651020N       | 4652020        |
| 20 µl    | 4651130N       | 4652130        |
| 25 µl    | 4651030N       | 4652030        |
| 50 µl    | 4651040N       | 4652040        |
| 100 µl   | 4651050N       | 4652050        |
| 200 µl   | 4651140N       | 4652140        |
| 250 µl   | 4651060N       | 4652060        |
| 500 µl   | 4651070N       | 4652070        |
| 1000 µl  | 4651080N       | 4652080        |
| 2000 µl  | 4651090N       | 4652090        |
| 3000 µl  | 4651100N       | 4652100        |
| 5000 µl  | 4651110N       | 4652110        |
| 10000 µl | 4651120N       | 4652120        |



Finnpipette F1 GLP Kit



Finnpipette F2 GLP Kit

## Finnpipette F1 and F2 Good Laboratory Pipetting (GLP) Kits

| Kit Number | Content                                                                                                                      | Finnpipette F1 | Finnpipette F2 |
|------------|------------------------------------------------------------------------------------------------------------------------------|----------------|----------------|
| GLP Kit 1  | • Single channel 1-10 µL 10-100 µL and 100-1000 µL pipettes<br>• F-Series Stand<br>• 3 Racks of Finntip Flex tips            | 4700850N       | 4700870        |
| GLP Kit 2  | • Single channel 0.2-2 µL, 2-20 µL, 20-200 µL and 100-1000 µL pipettes<br>• F-Series Stand<br>• 4 Racks of Finntip Flex tips | 4700860N       | 4700880        |
| GLP Kit 3  | • Single channel 10-100 µL, 100-1000 µL and 1-10 mL pipettes<br>• F-Series Stand<br>• 3 Racks of Finntip/Finntip Flex tips   | 4700865N       | 4700885        |
| GLP Kit 4  | • Single channel 2-20 µL, 20-200 µL and 100-1000 µL pipettes<br>• F-Series Stand<br>• 3 Racks of Finntip Flex                | 4701060N       | 4701070        |

## One family of pipettes is all you need to comfortably meet any application challenge

For decades, Thermo Scientific Finn timer pipetting systems have delivered unmatched productivity and ergonomics across thousands of labs and applications worldwide.

And the evolution of excellence continues with more enhancements – thanks to your valuable feedback – and it’s never been easier to select the perfect pipette to match your specific requirements.

### NEW Finn timer F1

The industry-leading manual pipette, thanks to its lightweight design, low plunger, tip attachment and ejection forces – plus a secure volume lock and adjustable finger rest.

### Finn timer F2

Rugged construction and ease of maintenance for your more demanding applications.

### Finn timer Novus

The move to electronic pipetting is easier than you think to support higher throughput.

### Finntip Pipette Tips

Engineered for maximum performance and comfort with every Finn timer.

