



2017

THERMOBLOCKS

www.falcinstruments.it



THERMOBLOCKS

- 5 | Heating Thermoblocks
- 5 | Cooling Thermoblocks

THERMOBLOCKS

Heating and cooling Thermoblocks



THERMOBLOCKS



TD: 25°C - 200°C



Ideal for storage, digesting, incubating, concentrating, and other heating processes. Useful for enzymatic and chemical-catalytic reactions. The temperature is set by three-step integrated programming

Heating Thermoblocks



TD 200 P1



TD 200 P2



TD 200 P2+



TD 200 P3

Technical data		TD 200 P1*	TD 200 P2*	TD 200 P2+*	TD 200 P3*
Blocks number	n°	1	2	2 separated	3
Adjustable temperature	°C	25 ÷ 200	25 ÷ 200	25 ÷ 200	25 ÷ 200
Operating temperature	°C	+5 room ÷ 200	+5 room ÷ 200	+5 room ÷ 200	+5 amb ÷ 200
Precision	°C	±0,5	±0,5	±0,5	±0,5
Timer	h:mm	9:59 + ∞	9:59 + ∞	9:59 + ∞	9:59 + ∞
Power	W	150	300	150+150	500
Dimensions LxWxH	mm	285x200x110	305x230x110	305x290x110	305x290x110
Weight	Kg	2,9	3,4	3,6	4,1
Code		621.0422.11	621.0422.12	621.0422.15	621.0422.13

V/HZ 220/240-50/60

CLASSE DI SICUREZZA 0

CLASSE DI PROTEZIONE IP 30

Cooling Thermoblocks



TDC: 0°C ÷ 100°C



Ideal for storage, digestion, incubation, concentration, and other heating processes. Useful for enzymatic reactions and chemical-catalytic. The temperature is set through three-step integrated programming



TDC 100 P1



TDC 100 P2

Technical data		TDC 100 P1	TDC 100 P2
Blocks number	n°	1	2
Adjustable temperature	°C	0 ÷ 100	0 ÷ 100
Operating temperature	°C	5 ÷ 100	5 ÷ 100
Precision	°C	±0,5	±0,5
Timer	h:mm	9:59 + ∞	9:59 + ∞
Program	step	3 temp + 3 time	3 temp + 3 time
Power	W	120	240
Dimensions LxWxH	mm	325x270x191	325x320x191
Weight	Kg	4,4	5,5
Code		621.0422.21	621.0422.22

V/HZ 220/240-50/60

CLASSE DI SICUREZZA 0


CLASSE DI PROTEZIONE IP 30

THERMOBLOCKS




One-place bucket with ceramic bead
dim 55x75x36. For TD models
Code 621.0423.50


Two place bucket with ceramic bead
dim 130x75x36. For TD models
Code 621.0423.51




6 holes - test tubes 28 mm.
For models TD
Code 621.0423.62




30 holes - test tubes 0.5 ml.
For TD models
Code 621.0423.72



96 holes - 0.2 ml - skirted PCR, vials,
double block version. For models TD
Code 621.0423.74




20 holes - test tubes
1.5 ml. For TD models
Code 621.0423.73




2 holes - Test tubes 50 ml flat bottom.
For TD models
Code 621.0423.64


THERMOBLOCKS




20 holes - test tubes 10 mm.
For TD models
Code 621.0423.56




72 holes - 0.2 ml - vials.
For models TD
Code 621.0423.70




20 holes - test tubes 13 mm.
For TD models
Code 621.0423.59



20 holes - test tubes 12 mm.
For TD models
Code 621.0423.58




12 holes - test tubes 16 mm.
For TD models
Code 621.0423.60




20 holes - test tubes 2 ml.
For models TD
Code 621.0423.57


THERMOBLOCKS




8 holes - test tubes 19 mm.
For TD models
Code 621.0423.61




4 holes - test tubes 33 mm.
For TD models
Code 621.0423.63




30 holes - test tubes 0.5 ml.
For TD models
Code 621.0423.55



96 holes - 0.2 ml - non skirted vials,
double block version.
For TD models
Code 621.0423.71




Elisa plate block.
For TD models
Code 621.0423.81




20 holes - test tubes 1.5 ml.
For models TDC
Code 621.0423.92


THERMOBLOCKS




30 holes - test tubes 0.5 ml.
For models TDC
Code 621.0423.93




96 holes - 0.2 ml - skirted PCR, vials,
double block version. For models TDC
Code 621.0423.95




20 holes - test tubes 1.5 mm.
For models TDC
Code 621.0423.91



20 holes - test tubes 10 mm.
For TDC models
Code 621.0423.85




72 holes - 0.2 ml - vials.
For models TDC
Code 621.0423.90




20 holes - test tubes 10 mm.
For TDC models
Code 621.0423.86

THERMOBLOCKS



96 holes - 0.2 ml - non skirted vials,
double block version. For models TDC
Code 621.0423.94



Elisa plate block.
double block version. For TDC models
Code 621.0423.98



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2 | 24047 Treviglio (BG) - Italy

+39 0363 304660 | falc@falcinstruments.it | falcinstruments.it