

LABNET INTERNATIONAL, INC.

MPS 1000™ Mini Plate Spinner

Finally, a compact and economical centrifuge for PCR plates

Labnet
Labnet International, Inc.



- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges

Patent Pending



SPECIFICATIONS

| | |
|---------------------|--|
| Capacity: | 2 PCR Plates |
| Rotor Speed: | 2500 rpm |
| G-Force: | 500 xg |
| Rotor: | Fixed Vertically |
| Dimensions (WxDxH): | 7.5 x 8.2 x 7.2 in. 19 x 21 x 18.3 cm |
| Weight: | 6 lbs / 3 kg |
| Electrical: | 120V or 230V, 50-60 Hz |

ORDERING INFORMATION

| Catalog No. | Description |
|-------------|--|
| C1000* | Mini Plate Spinner centrifuge with two position rotor |

* To order this product in 230V, please add (-230V) to the item number.

Labnet
Labnet International, Inc.

PO Box 841 • Woodbridge • NJ • 07095 • USA • website: <http://www.labnetlink.com>
US • toll free: 888-LABNET1 • fax: (732) 417-1750 • email: labnet@labnetlink.com
International: +1-732-417-0700 • fax: +1-732-417-1750 • international@labnetlink.com