

MILAB

Ultrasonic Baths Basic Line



Die casting stainless steel tank
Industrial grade integrated
circuit



20÷80°C temperature range



0÷20 minutes working time
Dewaxing and degreasing
function

MODEL	Tank Size (mm)	Overall Size (mm)	Capacity (L)	Frequency (KHz)	Ultrasonic Power (W)	Heating Power (W)	Timing (min)	Temp. (°C)	Drainage
UTA 1.3	150x140x65	190x170x195	1.3	40	60	100	0÷20/Always On	20÷80° C	None
UTA 3.2	240x140x100	270x170x240	3.2	40	120	100	0÷20/Always On	20÷80° C	None
UTA 6.5	300x155x150	330x180x330	6.5	40	180	300	0÷20/Always On	20÷80° C	Yes
UTA 10	300x240x150	330x270x330	10	40	240	300	0÷20/Always On	20÷80° C	Yes
UTA 15	330x300x150	360x330x330	15	40	360	400	0÷20/Always On	20÷80° C	Yes
UTA 22	500x300x150	500x330x330	22	40	480	500	0÷20/Always On	20÷80° C	Yes

FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2 | 24047 Treviglio (BG) - Italy | +39 0363 304660 | falc@falcinstruments.it | falcinstruments.it