

RUD

M1324

High-Speed Microcentrifuge

A Leading Manufacturer for Pre-clinical Research





M1324 High-Speed Microcentrifuge

Maximum capacity: 24 X 1.5/2.0 mL

Maximum speed: 21,130 x g (15,000 rpm)

Temperature control range: room temperature ,
RT mode and Cooling mode, one keyboard toggle

Automatic countdown display mode under unmanned status while spinning Highly intelligent analysis software: the LED display allows you to preprogram 12 programs simultaneously for directly use
Multiple abnormal alarm including over-speed, over-temperature, as well as unbalance, ensuring the safety of the operation of the user.

Rear-mounted compressor design, saving space

Technical parameters:

Temperature range	Room temperature
Max. speed	21,130 X g (15,000 rpm)
Max. capacity	24 X 1.5/2.0 mL
Display	4.3-inch touchscreen, 800*480 px
Dimensions (WXDXH)	300mm*514mm*297mm

RWD Life Science Co.,Ltd

Add: Room201, Robeta Building, No.11 High Tech North 4th Road, Nanshan District, Shenzhen City, Guangdong Province, 518057, China

Tel: (86)-755-86111286 ext 8283 Web: www.rwdstco.com E-mail: sales@rwdstco.com

RWD Life Science Inc.

Add: 6540 Lusk Blvd, Suite C161, San Diego, CA 92121

Tel: (858)900-5879 E-mail: sales@rwdstco.com