

# CertoClav Connect



## The fully automated including documentation

CertoClav Connect is a fully automated laboratory autoclave that requires no prior training. It is operated by means of a 7"- high resolution touch screen, the CertoTouch. All operating steps are explained with training videos directly on the touch screen. The software provides user management for comprehensive documentation of operation. Individual programs can be generated directly at the touch screen.

When the device is connected to the Internet via WLAN, it automatically connects to the CertoCloud via a secure line. This allows all protocols to be accessed via a browser, and customized programs to be created and sent to the devices. Programs created directly on the device are synchronized as well.

Thanks to the automated documentation of users and cycles, the CertoClav Connect meets all GLP requirements.

## Important features

- fully automated
- automatic venting
- electric exhaust valve
- optional flexible media sensor
- Time and temperature can be selected individually
- Controller with LCD display
- Venting by fractional flow procedure
- Ethernet, LAN, Wi-Fi
- 7" tablet with video assistant
- Cloud Connection
- Heating integrated in the base
- electronic safety interlock

### Scope of delivery:

incl. instrument tray,  
power cord, 7" Tablet



**PT100-Temperatursensor –**  
For autoclaving liquids  
optionally available.



## Specifications

	CertoClav Connect
Article number	8510176
Chamber volume	18 l
Chamber diameter	26 cm
Chamber depth	35 cm
Chamber material	stainless steel 316L
Max. operating pressure	2,7 bar
Mains connection	230 V, 50/60 Hz
Heating capacity	1300 W
Max. sterilization temperature	20 °C bis 140 °C, einstellbar
Net weight	16 kg
Product dimensions (L*W*H)	44*38*54 cm

## Accessories

Article	Article number
Wire basket with handle	1300149
Service Set 2	1310186
Media sensor	1300185
CertoServer	8500508
CertoCloud premium license	8500212
Distiller for water	8500482

Factory calibration possible.