

RBS 50



Liquid cleaning agent, concentrated, alkaline and non-foaming. RBS 50 is suitable for use in ultrasonic baths and automatic washing machines.

Advantages

- Cost effective
- Efficient at low concentrations
- High cleaning and degreasing power
- Multipurpose detergent
- Alternative to solvents and sulfochromic acids known as toxic and hazardous
- Contains chlorine that provide decontaminating properties
- Total rinse



Use

In laboratories, RBS 50 is used for cleaning glassware and laboratory equipment made of stainless steel, plastic, porcelain, etc. In industry, it is recommended for cleaning parts and surfaces made of glass, stainless steel, production equipment, etc.

Properties

- pH pure product : 13.5
- pH 0.3% v/v in water : 10.4

Concentration

Use concentration : 0.3 to 0.5% v/v in water.

Packing

RBS 50 is available in 5 l blue polyethylen bottle (only available in box of 4 x 5 l).