



WASHING SOLUTIONS
for pharma sector



PRECISE RESULTS AND GUARANTEED COSTS

Smeg SpA, the professional washing market leader, offers a range of solutions for the correct treatment of pharmaceutical punches or tools for production departments. The only answer to the need for certain, safe and repeatable results.

PERFECT WASHING

- Specific trolley to wash dies.
- custom-made programs for both an accurate yet delicate wash.
- specific solutions for different needs



OPTIMAL DRYING

- High efficiency drying system.
- Integrated hepa filter for clean uncontaminated air.
- The goods must not be manually washed.

ERGONOMICS

Operations performed safely, avoiding improper lifting of weights by staff



TIME SAVING

Optimization of processing time, and management of human resources

CUSTOM-MADE SERVICES

- IQ/OQ/PQ on demand
- custom-made after-guarantee contracts. Servicing and repairs by specialized approved technicians.
- tailor-made solutions for the specific requirements of the customers.



SECURE INVESTMENT

- 24 months warranty
- impeccable results
- high profitability of the production process
- saving time and reducing consumption
- optimized working processes

PW RANGE

60cm or 90cm wide washing machines with an ergonomic and safe base.



- External dimensions 600x605x850 mm
- Tank capacity: 150lt
- Free-standing positioning
- Possibility to install under worktop
- Printer
- Three-phase power supply
- Custom-made and specific wash programs for each type of utensil
- Ergonomic base



- External dimensions 900x605x850 mm
- Tank capacity: 150lt
- Free-standing positioning
- Possibility to install under worktop
- Detergent compartment:
- Integrated printer
- Three-phase power supply
- Custom-made and specific wash programs for each type of utensil
- Ergonomic base

PW RANGE

Large washing machines for vast loads and volumes of work



- External dimensions 650x720x1945 mm
- Tank capacity: 266lt
- Drop down door in glass
- Detergent compartment:
- Integrated printer
- Three-phase power supply
- Custom-made and specific wash programs to each type of utensil




- External dimensions 780x800x1945 mm
- Tank capacity: 266lt
- Sliding glass door
- Detergent compartment:
- Integrated printer
- Three-phase power supply
- Custom-made and specific wash programs for each type of utensil

WORKING TOGETHER TO DELIVER THE BEST WASHING SOLUTIONS

Designed and built to measure to respond to differing washing needs.

Documentation and testing according to user requirement



Process traceability and precise results

Process optimization: human resources efficiently used for high productivity.



Smeg S.p.A.
Via Leonardo da Vinci, 4
42016, Guastalla (RE), Italia
instruments@smeg.it
www.smeg-instruments.com