

Elmasonic Accessory - Inserts

Cleaning glasses with plastic cover and rubber ring



- can be used for simultaneous cleaning of parts with different chemicals, for removing abrasive substances or for cleaning with aggressive chemicals (e.g. acids or alkalis)
- can either be placed directly in a cleaning basket or in combination with an insert cover be hung in the ultrasonic bath




	600 ml	1000 ml
Order number	1047272	1047376
Glas with cover dimensions ext. Dia/H/D (mm)	110/135/125	110/185/125
Glas dimensions internal Dia _i /H _i (mm)	80/105	85/160
Rubber ring external diameter Dia _R (mm)	100	

Dimension definition

Cleaning glas
with cover



Rubber ring



All prices ex works, excluding packaging, plus VAT.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.