


Printing & Viewing Tips

INDIVIDUAL CANISTERS, FRAMES AND PLATFORM DIVIDERS

Canisters constructed of lightweight aluminum including locking handle and thumb hole. Frames constructed of corrosion resistant stainless steel including lifting handle. Custom sizes available upon request.

INDIVIDUAL CANISTERS



BAXTER FENWAL

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
4R9951	ZC060H	3 ⁵ / ₈ " x 6 ¹ / ₄ " x 3 ⁵ / ₈ "
4R9953	ZC021	6" x 7 ¹ / ₂ " x 3 ⁵ / ₈ "
4R9955	ZC022	6" x 9 ⁵ / ₈ " x 3 ⁵ / ₈ "
4R9957	ZC025	6" x 11" x 3 ⁵ / ₈ "
4R9959	ZC023	6" x 14" x 3 ⁵ / ₈ "
4R2111	ZC024	5 ³ / ₄ " x 8" x 5 ⁵ / ₈ "


CHARTERMED

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
7005-2	ZC004	7 ¹ / ₄ " x 9" x 3 ⁵ / ₈ "
7039-2	ZC004	7 ¹ / ₄ " x 9" x 3 ⁵ / ₈ "
7303-2	ZC004	7 ¹ / ₄ " x 9" x 3 ⁵ / ₈ "
7400-2	ZC005	9" x 9" x 3 ⁵ / ₈ "
7410-2	ZC005	9" x 9" x 3 ⁵ / ₈ "
7403-2	ZC005	9" x 9" x 3 ⁵ / ₈ "
SPECIAL	ZC002	11 ¹ / ₂ " x 11" x 3 ⁵ / ₈ "
PEDIATRIC	ZC001	5 ³ / ₄ " x 9 ¹ / ₈ " x 7 ¹ / ₁₆ "

GAMBRO

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
DF-200	ZC011	8" x 6 ³ / ₁₆ " x 1 ³ / ₁₆ "
DF-700	ZC012	12" x 6 ³ / ₁₆ " x 1 ³ / ₁₆ "
DF-1000	ZC013	14 ¹ / ₄ " x 6 ³ / ₁₆ " x 1 ³ / ₁₆ "
DF-1200	ZC014	16 ¹ / ₂ " x 6 ³ / ₁₆ " x 1 ³ / ₁₆ "


INDIVIDUAL FRAMES



<u>FRAME</u>	<u>PART NUMBER</u>
2 PLACE	ZF2
3 PLACE	ZF3
4 PLACE	ZF4
5 PLACE	ZF5
6 PLACE	ZF6
7 PLACE	ZF7
8 PLACE	ZF8

NOTE: To order, add canister part number after frame part number. Example: ZF4-ZC021

INDIVIDUAL PLATFORM DIVIDERS



<u>FREEZER DIAMETER</u>	<u>PART NUMBER</u>
18 INCH	ZD18
19 INCH	ZD19
20 INCH	ZD20
31 INCH	ZD31
39 INCH	ZD39

NOTE: To order, add canister part number after platform divider part number. Example: ZD31-ZC021

Printing & Viewing Tips
COMPLETE SYSTEMS
CANISTERS, FRAMES AND PLATFORM DIVIDERS

Systems consist of complete set of canisters, frames plus platform divider. Canisters constructed of lightweight aluminum including locking handle and thumb hole. Frames including lifting handle and platform dividers are constructed of corrosion resistant stainless steel. Custom sizes available upon request.

20 INCH FREEZER

BAXTER FENWAL

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	434	62 (7 PL)	ZS200	27.75"
4R9953	224	56 (4 PL)	ZS201	26"
4R9955	164	41 (4 PL)	ZS202	26"
4R9957	136	34 (4 PL)	ZS203	26"
4R9959	88	22 (4 PL)	ZS204	26"

CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	129	43 (3 PL)	ZS205	24"
7039-2	129	43 (3 PL)	ZS205	24"
7303-2	129	43 (3 PL)	ZS205	24"
7400-2	129	43 (3 PL)	ZS206	27.75"
7410-2	129	43 (3 PL)	ZS206	27.75"
7403-2	129	43 (3 PL)	ZS206	27.75"
SPECIAL	70	35 (2 PL)	ZS207	25"

GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	108	27 (4 PL)	ZS208	27"
DF-700	64	16 (4 PL)	ZS209	27"

39 INCH FREEZER

BAXTER FENWAL

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	1932	276 (7 PL)	ZS900	27.75"
4R9953	912	228 (4 PL)	ZS901	26"
4R9955	752	188 (4 PL)	ZS902	26"
4R9957	608	152 (4 PL)	ZS903	26"
4R9959	456	114 (4 PL)	ZS904	26"

CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	552	184 (3 PL)	ZS905	24"
7039-2	552	184 (3 PL)	ZS905	24"
7303-2	552	184 (3 PL)	ZS905	24"
7400-2	552	184 (3 PL)	ZS906	27.75"
7410-2	552	184 (3 PL)	ZS906	27.75"
7403-2	552	184 (3 PL)	ZS906	27.75"
SPECIAL	304	152 (2 PL)	ZS907	25"

GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	464	116 (4 PL)	ZS908	27"
DF-700	320	80 (4 PL)	ZS909	27"

31 INCH FREEZER

BAXTER FENWAL

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	1120	160 (7 PL)	ZS100	27.75"
4R9953	560	140 (4 PL)	ZS101	26"
4R9955	440	110 (4 PL)	ZS102	26"
4R9957	352	88 (4 PL)	ZS103	26"
4R9959	256	64 (4 PL)	ZS104	26"

CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	336	112 (3 PL)	ZS105	24"
7039-2	336	112 (3 PL)	ZS105	24"
7303-2	336	112 (3 PL)	ZS105	24"
7400-2	336	112 (3 PL)	ZS106	27.75"
7410-2	336	112 (3 PL)	ZS106	27.75"
7403-2	336	112 (3 PL)	ZS106	27.75"
SPECIAL	180	90 (3 PL)	ZS107	27"

GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	288	72 (4 PL)	ZS108	27"
DF-700	192	48 (4 PL)	ZS109	27"

