



LAB 500 Series

Undercounter Glassware Washers




- Chamber volume ~171 lt / 6.04 ft³
- Basket volume ~151 lt / 5.33 ft³



LAB 500 SC

- Washing system **on two independent levels**
- Chamber drying by electrical heating elements
- Stainless steel or high visibility full glass door



LAB 500 CL

- Washing system **on two independent levels**
- Forced hot air drying system **on two independent levels**



Stainless steel door version **LAB 500 SC** only

Control panel

- LED display, 10 programs for laboratory glassware: 5 pre-programmed cycles, 5 customizable cycles.

Connections

- RS 232 port dedicated for printer or PC connection to monitor and validate the washing cycles and/or the data storage.



Stainless steel door version **LAB 500 SC** and **LAB 500 CL**

Control panel

- LCD display, 40 programs for laboratory glassware: 20 pre-programmed cycles, 20 customizable cycles.

Connections

- RS 232 port for printer or PC connection.
- USB port on front panel for cycle data storage and the program updating.



Full glass door version **LAB 500 SC**, **LAB 500 CL** and **LAB 500 DRS**

Control panel

- Full glass, soft touch control panel, LCD display, 40 programs for laboratory glassware: 20 pre-programmed cycles, 20 customizable cycles.

Connections

- RS 232 port for printer or PC connection.
- USB port on front panel for cycle data storage and the program updating.